

Med-Aire Alternating Pressure Mattress Replacement System with Low Air Loss

14027

Mattress and Pump System,
80" (L) x 36" (W) x 8" (H), 1/ea

HCPCS code: E0277

FEATURES:

- This combination therapy system provides both alternating pressure and low air loss to optimize pressure redistribution and manage skin maceration.
- Designed to prevent, treat and manage pressure ulcers in the home, long term care, or acute care setting.
- **2-way stretch cover** is low shear, fluid-resistant, vapor permeable, quilted and zippered.
- Each one of the 20, **8" deep air bladders** are easily removed and replaced.
- Pillow function maintains the 3 air cells at the head of the bed in static mode for patient comfort.
- 3 diagnostic alarms tell you whether a problem is with the mattress or the pump for more efficient servicing.
- 350 lb. weight capacity.
- CPR valve allows for rapid deflation.

PARTS:

- 14027M Mattress Only, 1/ea
- 14027P Pump Only, 1/ea
- 14027-CELL Replacement Bladder for 14027 and 14028N, 1/ea
- 14027-CELL 1 Replacement Bladder with Laser Holes for 14027 and 14028N, 1/ea

SPECIFICATIONS:

Pump:

- Power 120 VAC, 60 HZ..
- Pressure range 30 mmHg to 80 mmHg (inner mattress pressure).
- Air volume 8 LPM.
- Cycle Time: 10 minutes.
- Functions: Static and Dynamic Features.
- Power Cord Length: 15 Ft.
- Pump Dimensions: 11" (W) x 8" (H) x 3.5" (D).
- UL: 60601-1, 1st Edition-2006.
- Limited 1 Year Warranty.

Available
1 Pack

OUR MOST COST-EFFECTIVE PRESSURE REDISTRIBUTING GROUP II MATTRESS REPLACEMENT!

MATTRESS CAN OPERATE IN AN ALTERNATING OR STATIC PRESSURE MODE!

350 LB
WEIGHT CAPACITY



14427 (Electric Pump)

14027

Mattress:

- Dimensions: 80" (L) x 36" (W) x 8" (H).
- Base Material: zippered PVC mesh.
- Maximum Support Capacity: 350 lbs.
- Limited 1 Year Warranty.



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

PRESSURE PREVENTION

drive

To order call toll free: 877.224.0946

13