

COMPARE TO:
drive™ **EHOB®**
14428 **1006EC**



14428



14427 (Electric Pump)



14426 (Manual Pump)

Available 1 Pack

Premium Guard - Static Air Mattress Overlay

- 14428 72" (L) x 35" (W) x 2.5" (H), Static Air Mattress, 1/ea
- 14426 Manual Pump for Premium Guard Air Mattress, 1/ea
- 14427 Electric Pump for Premium Guard Air Mattress, 1/ea

HCPCS code: E0197

FEATURES:

- Unique construction allows for efficient air flow throughout mattress for use in the treatment and prevention of pressure sores.
- Easy to clean with mild detergent.
- Emergency air release valve deflates mattress quickly and safely.
- Air flow holes reduce heat and fluid build up.
- Adjustable bed straps prevent mattress from moving.
- Two adapters for easy filling and air removal.
- 250 lb. weight capacity.
- Limited 90 Day Warranty on seams and valve.



14400



NEW PACKAGING!

Fill Hose



14401

Available 1 Pack

Premium Guard Water Mattress

- 14400 68" (L) x 32" (W) x 4" (H), Water Mattress, 1/ea
- 14401 Fill hose for Water Mattress (Sold separately)

FEATURES:

- Three water baffles enhance movement and support to reduce pressure and conform to varying body weights.
- Filling lines and "foot placement" lines provide for proper patient positioning.
- Double sealed seams provide maximum strength and durability.
- Vinyl construction is durable and easy to clean with mild detergent.
- Fill lines are conveniently located.
- 250 lb. weight capacity.
- Limited 90 Day Warranty on seams and valve.

SPECIFICATIONS

Item #	Length	Width	Height	Product Weight	Carton Shipping Weight	Weight Capacity
14428	72"	35"	2.5"	8 lbs.	11 lbs.	250 lbs.
14400	68"	32"	4"	4 lbs.	6 lbs.	250 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.