

Folding, Portable, Upholstered Commode with Wheels and Drop Arm

Available
1 Pack


CHROME MODEL SHOWN

COMES STANDARD WITH PADDED SEAT.

TOOL FREE ASSEMBLY



ALSO AVAILABLE IN SILVER VEIN FINISH



11120SV-1 SHOWN WITH FOOTRESTS (NOT INCLUDED)




- 11120SV-1 Assembled, Silver Vein, (Has replaced 11120-1), 1/cs
- 11120KD-1 Knock Down, Tool Free Assembly, Chrome 1/cs
- PSTDSF-TF Swing-away Tool Free Footrests, (Plastic) (For use on 11120SV-1), 1 pr/cs
- 11120B Commode Bucket for 11120, 1/cs
- 11120-SP Replacement Seat, 1/cs
- 11120-S Replacement Seat Pad, 1/cs
- STDSF-TF Swing-away Footrests Tool Free, (Metal) (For use on 11120KD-1), 1 pr/cs

FEATURES:

- Chrome plated or Silver Vein coated steel tubing frame with welded construction.
- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.
- Swing-away and detachable footrests (sold separately).
- Injection molded large plastic toilet seat, 12 quart bucket and cover come standard.
- Removable, tool free back.
- Rust resistant, swivel casters (2 with locks).
- Can be shipped via UPS.
- Limited Lifetime Warranty.

11120KD-1



NEW REINFORCED TOOL FREE BACK REST DESIGN



3" NON-SKID, RUST RESISTANT, SWIVEL CASTERS. TWO WITH LOCKS!



Available
1 Pack

K.D. Aluminum Shower Chair/Commode with Casters

- 11114KD-1 Knock Down, (Has replaced 11114-2), 1/cs
- 11114SEAT Replacement Seat, 1/cs

FEATURES:

- Easy to assemble frame.
- Padded, open front, vinyl toilet seat.
- Aluminum construction is lightweight and durable.
- Non-skid, rust resistant, swivel casters (2 rear locking).
- Removable, tool free back.
- Plastic arms for added comfort.
- Comes with removable 12 qt. commode bucket and cover.
- Limited Lifetime Warranty.

NEW IMPROVED CASTERS!



See Page 15
Biodegradable Sanitary Bag Commode Liner
RTL12085
12 Commode Liners Per Box, 6bx/cs

EASILY ASSEMBLED WITHOUT TOOLS!



SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11120KD-1	18"	22"	17"	17"	20"	20"	35 lbs.	39 lbs.	350 lbs.
11120SV-1	18"	22"	17"	17"	20"	20"	35 lbs.	39 lbs.	350 lbs.
11114KD-1	18" (8" ARM HEIGHT)	21.5"	16"	16"	21"	21"	14.10 lbs.	18.1 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.