



RS185001

RS185003

RS185005

FRONT CASTERS COME WITH "PIN LOCKS" THAT CAN BE SET IN THE SWIVEL OR LOCKED POSITION!



COMPARE TO:	
drive	Graham-Field®
RS185001	12022010
RS185005	12022020
RS185002	12022011
RS185003	12022013



RS185004S



11124

SPECIFICATIONS

Item #	Overall Width Open	Overall Width Closed	Overall Height	Seat To Floor Height	Back Height	Width Seat	Depth Seat	Overall Length (with riggings)	Product Weight	Carton Shipping Weight	Weight Limit
RS185001	24.5"	11"	37.5"	22.5"	15"	17"	18"	42"	51 lbs.	56 lbs.	300 lbs.
RS185002	24.5"	11"	37.5"	22.5"	15"	17"	18"	42"	51 lbs.	56 lbs.	300 lbs.
RS185003	23"	9"	37.5"	22.5"	15"	17"	18"	34.75"	48 lbs.	53 lbs.	300 lbs.
RS185005	24.5"	11"	37.5"	22.5"	15"	17"	18"	42"	56 lbs.	63 lbs.	300 lbs.

Stainless Steel Rehab/Shower Commode

Available
1 Pack

- RS185001 With Removable Full Arms, 24" Wheels and Swing-away Footrests, 1/cs
- RS185002* With Removable Full Arms, Swing-away Footrests, Four 5" Casters, Front Caster Pinlocks, 1/cs
- RS185003 With Removable Full Arms, Swing-away Footrests, Four 5" Casters, Rear Casters with Locks 1/cs
- RS185005 With Removable Full Arms, 24" Wheels, Swing-away Footrests and Open Seat, 1/cs
- PSTDSF-TF Replacement Swing-away Footrest, Tool Free (Plastic footplate) (All rehab/shower commodes come with riggings), 1/cs
- RS185003S Replacement Seat (closed) for RS185001/RS185003 and RS185002, 1/ea
- RS185004S Replacement Seat (open) for RS185005, 1/ea
- 11124 Replacement Commode Bucket for RS Series Rehab/Shower Commode, 1/ea

FEATURES:

- Can be positioned over a standard toilet or used as a portable self contained commode.
- Removable full arms with plastic armrests allow for easy transfer.
- Corrosion resistant stainless steel frame.
- Cushioned, vinyl seat is sealed to prevent moisture penetration.
- Removable pail and lid are included.
- Mag-style 24" rear wheel with rubber tires allows individual to self propel. (RS185001, RS185005).
- RS185003 features four 5" casters, rear casters with locks.
- Swing-away, detachable footrests with durable, non-corrosive footplates come standard.
- Warranty: Lifetime on Frame.



See Page 15
Biodegradable Sanitary Bag Commode Liner
RTL12085
12 Commode Liners Per Box, 6bx/cs

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.