

18" WIDTH BETWEEN FRONT ARMS!



RTL12027RA



REMOVABLE ARMS ALLOW YOU TO STOCK ONE PRODUCT AND FILL THE NEEDS OF TWO.

REMOVEABLE ARMS PINLOCK INTO METAL SLEEVES!

18" WIDTH BETWEEN ARMS!



12008KDR

TOOL FREE REMOVABLE ARMS!



2 in 1 Locking Elevated Toilet Seat with Tool Free Removable Arms

Available **1 Pack**

RTL12027RA

With Removable Arms, Retail Packaged, (Has replaced 12027RA), 1/cs

12027RA-4BULK

With Removable Arms, 4/cs (Each toilet seat individually polybagged, 4 per carton)

FEATURES:

- Designed for individuals who have difficulty sitting down, or getting up from the toilet.
- Tool Free Arms can be removed or added as needed.
- Arm sleeves and receivers are made of metal for more durability.
- Heavy-duty molded plastic construction provides additional strength and durability.
- Improved locking device is engineered with larger, heavy-duty "worm screw" and locking plate to ensure a safe, secure fit.
- Lightweight and portable.
- Easily attaches to toilet bowl.
- No tools required for installation.
- Easy to clean.
- Fits most toilets.
- Limited Lifetime Warranty.



K.D. 2 in 1 Locking, Elevated Toilet Seat with Tool Free Removable Arms

Available **1 Pack**

12008KDR

(Has replaced 12008), 1/cs

FEATURES:

- Designed for individuals who have difficulty sitting down, or getting up from the toilet.
- Tool Free Arms can be removed or added as needed.
- Arm sleeves and receivers are made of metal for more durability
- Heavy-duty molded plastic construction provides additional strength and durability.
- Standard locking mechanism ensures safety.
- Lightweight and portable.
- Easily attaches to toilet bowl.
- No tools required for installation.
- Easy to clean.
- Fits most toilets.
- Limited Lifetime Warranty.



SPECIFICATIONS

Item #	Width (seat)	Depth (seat)	Seat Height	Width Between Arms	Product Weight	Carton Shipping Weight	Weight Limit
12027RA	17"	16.5"	5"	18"	5 lbs.	8.5 lbs.	300 lbs.
12027RA-4BULK	17"	16.5"	5"	18"	5 lbs.	23 lbs.	300 lbs.
12008KDR	15.5"	16"	5"	18"	4.5 lbs.	6.5 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.