



12440-2

**ENVIRONMENT FRIENDLY PRODUCT**

Recycled materials were used  
to make this product.

### Bath Bench with Removable Padded Arms

Available  
**1 Pack**

12440-2 2/cs

12440-1 1/cs

12002T Replacement Tips for 12440, 12445, 4/bx

**FEATURES:**

- Removable padded armrests provide comfort and added security.
- Blow-molded bench provides comfort and strength.
- Drainage holes in bench reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Adjustable height legs.
- Tool free assembly.
- Limited Lifetime Warranty.



12445-2

**ENVIRONMENT FRIENDLY PRODUCT**

Recycled materials were used  
to make this product.

### Bath Bench with Back and Removable Padded Arms

Available  
**1 Pack**

12445-2 (Has replaced 12202A-4), 2/cs

12445-1 1/cs

12002T Replacement Tips for 12440, 12445, 4/bx

**FEATURES:**

- Removable padded armrests provide comfort and added security.
- Blow-molded bench and back provide comfort and strength.
- Drainage holes in bench reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Adjustable height legs.
- Tool free assembly.
- Limited Lifetime Warranty.

**SPECIFICATIONS**

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat (adjustable)	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12440-2	16"	16"	17.3"	14" (Fixed)	20"	16"	6.75 lbs.	15 lbs.	300 lbs.
12445-1/2	16"	16"	17.3"	14", 15", 16.5"	20"	16"	9.68/19.1 lbs.	10.1/21.3 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.