

SPECIFICATIONS

| Item # | Width (outside legs) | Width Seat | Depth (outside legs) | Depth Seat | Seat Height (highest) | Seat Height (lowest) | Product Weight | Carton Shipping Weight | Weight Limit |
|-------------------|----------------------------|---------------|----------------------------|---------------|-----------------------------|----------------------------|-------------------|------------------------------|-----------------|
| 12105KDR-2 | 19" | 19.75" | 18" | 12" | 19.75" | 15.75" | 8.2 lbs. | 19 lbs. | 300 lbs. |
| 12105-4 | 19" | 19.75" | 18" | 12" | 19.75" | 15.75" | 8.2 lbs. | 37 lbs. | 300 lbs. |