



RTL12105KDR

NEW PINCH FREE COVER ALLOWS FOR PUSH PINS TO BE DEPRESSED WITHOUT PINCHING FINGER. THE AUDIBLE "CLICK" WILL BE HEARD WHEN SECURE!

RETAIL MODEL COMES WITH CONVENIENT CARRY BAG.



Replacement Carry Bag 12105B



ENVIRONMENT FRIENDLY PRODUCT
Recycled materials were used to make this product.



**READY SET GO
K.D. Deluxe, Aluminum,
Bath Bench with Carry Bag**

RTL12105KDR With Back, Knock Down, Retail, Comes with Carry Bag, (Has replaced 12105KDR-2), 2/cs Replacement Tips for 12105, 4/bx
12002T Replacement Ready Set go Bag, 1/ea
12105B

FEATURES:

- Easy, tool free assembly of legs and back.
- Blow-molded bench and back provides comfort and strength.
- Drainage holes in seat and back reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Height adjustable, angled legs with suction style tips provide additional stability.
- Support collar prevents leg movement.
- Comes with a durable, nylon carry bag for easy and convenient transport.
- Limited Lifetime Warranty.



**ASSEMBLED WITH TOOL FREE
REMOVABLE BACK**

LEGS AND BACK EASILY ASSEMBLED WITHOUT TOOLS!



SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
RTL12105KDR	19"	19.75"	18"	12"	18"	14"	8 lbs.	19 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.