

SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12203KD-4	17"	20"	17"	12"	19.75"	15.75"	4.5 lbs.	21 lbs.	400 lbs.
12203KDR-1/2	17"	20"	17"	12"	19.75"	15.75"	5 lbs.	7/12 lbs.	400 lbs.
12203KDRB, T, Y-1	17"	20"	17"	12"	19.75"	15.75"	5 lbs.	7 lbs.	400 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.