



SUPPORT COLLAR PREVENTS LEG MOVEMENT!
12203KD-4

NEW PINCH FREE COVER ALLOWS FOR PUSH PINS TO BE DEPRESSED WITHOUT PINCHING FINGER. THE AUDIBLE "CLICK" WILL BE HEARD WHEN SECURE!



ENVIRONMENT FRIENDLY PRODUCT
Recycled materials were used to make this product.

Deluxe, K.D. Aluminum, Bath Bench without Back **Available 1 Pack**

- 12203KD-4 Knocked Down, Tool Free Assembly, (Has replaced 12106-4 and 12203-4), 4/cs
- RTL12203KDR Knocked Down, Tool Free Assembly, Retail, (Has replaced 12203KDR-1), 1/cs
- 12203KDRB-1 Knocked Down, Tool Free Assembly, Retail, Blue, 1/cs
- 12203KDRT-1 Knocked Down, Tool Free Assembly, Retail, Teal, 1/cs
- 12203KDRY-1 Knocked Down, Tool Free Assembly, Retail, Yellow, 1/cs
- 12002T Replacement Tips for 12203, 4/bx

FEATURES:

- Easy snap together assembly, Seat and Legs.
- Blow-molded bench provides comfort and strength.
- Drainage holes in bench reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Height adjustable, angled legs with suction style tips provide additional stability.
- Support collar reduces rattle.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12203KD-4	17"	20"	17"	12"	19.75"	15.75"	4.5 lbs.	21 lbs.	400 lbs.
RTL12203KDR	17"	20"	17"	12"	19.75"	15.75"	5 lbs.	7 lbs.	400 lbs.
12203KDRB, T, Y-1	17"	20"	17"	12"	19.75"	15.75"	5 lbs.	7 lbs.	400 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.