

HANDLE WIDTH IS ADJUSTABLE!



RTL12087



RTL12079



**TOOL FREE ASSEMBLY!**

**COMPARE TO:**  
**drive™ CAREX® MOEN®**  
**RTL12087 B368 DN7015**

**Toilet Safety Rail** **Available 1 Pack**  
 RTL12087 2 retail boxes/cs



**FEATURES:**

- Allows individuals to safely sit or stand.
- Easy assembly and installation.
- Comes with soft, padded handles for comfort.
- Width adjusts from 16" - 18" between handles.
- Limited Lifetime Warranty.



**Stand Alone Toilet Safety Rail** **Available 1 Pack**  
 RTL12079 2 retail boxes/cs

**FEATURES:**

- Free standing around toilet.
- Easy, tool free assembly using wing nuts.
- Fits standard and elongated toilets.
- Comes with padded armrests for added comfort.
- Comes with a magazine rack for added convenience.
- Limited Lifetime Warranty.

**SPECIFICATIONS**

Item #	Width Between Arms (smallest)	Width Between Arms (largest)	Depth	Height Seat	Product Weight	Carton Shipping Weight	Weight Limit
RTL12087	19.5"	22.5"	16"	11.5"	10.5 lbs.	24 lbs.	300 lbs.
RTL12079	NA	26.5"	19.5"	NA	9.0 lbs.	25.5 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.