



Available
1 Pack

K.D. Padded Transfer Bench

- 12005KD-2 Knock Down Padded, Tool Free Assembly Back, Legs and Arm, (Has replaced 12005-2), 2/cs
- 12005KD-1 Knock Down Padded, Tool Free Assembly Back, Legs and Arm, (Has replaced 12005-1), 1/cs
- 12005KDR-1 Knock Down Padded, Tool Free Assembly Back, Legs and Arm, Retail, 1/cs
- 12011SS Replacement Suction Cup for Transfer Bench, Small, 1pr/bx
- 12011SL Replacement Suction Cup with Screws for Transfer Bench, Large, 1pr/bx
- 12005TP Replacement Transfer Pad, 1/ea

FEATURES:

- Comfortable cushioned seat and backrest.
- "A" frame construction provides additional stability.
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs.
- All transfer benches are reversible to accommodate any bathroom.
- Aluminum frame is lightweight, sturdy and corrosion resistant.
- No exposed hardware to injure patient.
- Extra large suction cups provide added safety.
- Tool-Free Assembly requires connecting back, legs and arm.
- Limited Lifetime Warranty.


SUPPORT COLLAR PREVENTS LEG MOVEMENT!




MUST PULL "TAB" ON SUCTION CUP TO RELEASE FROM THE TUB.

12005KD-2

NEW PINCH FREE COVER ALLOWS FOR PUSH PINS TO BE DEPRESSED WITHOUT PINCHING FINGER. THE AUDIBLE "CLICK" WILL BE HEARD WHEN SECURE!



ASSEMBLED WITH TOOL FREE REMOVABLE BACK



SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12005KD-2	27.5"	24"	16.5"	16"	21.75"	17.75"	11.75 lbs.	26.5 lbs.	400 lbs.
12005KD-1	27.5"	24"	16.5"	16"	21.75"	17.75"	11.75 lbs.	14 lbs.	400 lbs.
12005KDR-1	27.5"	24"	16.5"	16"	21.75"	17.75"	11.75 lbs.	16 lbs.	400 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.