



TOOL FREE ASSEMBLY!

TOOL FREE REVERSABLE BACK!



MUST PULL "TAB" ON SUCTION CUP TO RELEASE FROM THE TUB.
RTL12011KDR

SUPPORT COLLAR PREVENTS LEG MOVEMENT!

NEW PINCH FREE COVER ALLOWS FOR PUSH PINS TO BE DEPRESSED WITHOUT PINCHING FINGER. THE AUDIBLE "CLICK" WILL BE HEARD WHEN SECURE!

12011KDR-1 GREAT FOR A PLANOGRAM!



ENVIRONMENT FRIENDLY PRODUCT
Recycled materials were used to make this product.

Available 1 Pack

K.D. Plastic Transfer Bench

- 12011KD-2 Knock Down Tool Free Back, Legs and Arm, 2/cs (Has replaced 12011-2 and 12020-2)
- 12011KD-1* Knock Down Tool Free Back, Legs and Arm, 1/cs (Has replaced 12011-1)
- RTL12011KDR Knock Down, Retail, Tool Free Back, Legs and Arm, (Has replaced 12011KDR-1), 1/cs
- 12011SS Replacement Suction Cup for Transfer Bench, Small, 1pr/bx
- 12011SL Replacement Suction Cup with Screws for Transfer Bench, Large, 1pr/bx
- 12011-A Replacement Arm for 12011KD and 12011 1/ea

FEATURES:

- "A" frame construction provides additional stability.
- Durable blow-molded plastic bench and backrest.
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs.
- Aluminum frame is lightweight, sturdy and corrosion resistant.
- All transfer benches are reversible to accommodate any bathroom.
- No exposed hardware to injure patient.
- Extra large, locking suction cups provide added safety.
- Tool-Free Assembly requires connecting back, legs and arm.
- Limited Lifetime Warranty.

LEGS AND BACK EASILY ASSEMBLED WITHOUT TOOLS!



SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12011KD-2	33"	26"	24"	18.5" - 19.5"	21.75"	17.75"	10.5 lbs.	24 lbs.	400 lbs.
12011KD-1	33"	26"	24"	18.5" - 19.5"	21.75"	17.75"	10.5 lbs.	13 lbs.	400 lbs.
RTL12011KDR	33"	26"	24"	18.5" - 19.5"	21.75"	17.75"	10.5 lbs.	14 lbs.	400 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.