



Portable non-permanent track system

X-Y Semi-permanent rack

Our KWIKtrak X-Y semi-permanent rack is a robust free standing room covering system with the capacity to handle everything including your bariatric lifting tasks

Suitable for use with fixed ceiling lifts or one of our portable lifts. Our semi-permanent rack is the ideal solution for temporary yet portable heavy duty lifting.


It is installed with a 180 mm KWIKtrak and can offer a maximum lifting capacity of 600 lb (270 kg).



Description			Description			Description		
# Product	Measure 1	Measure 2	# Product	Measure 1	Measure 2	# Product	Measure 1	Measure 2
9903030	3.0 m	3.0 m	9903540	4.0 m	3.5 m	9904045	4.5 m	4.0 m
9903035	3.5 m	3.0 m	9903545	4.5 m	3.5 m	9904545	4.5 m	4.5 m
9903040	4.0 m	3.0 m	9904040	4.0 m	4.0 m			
9903045	4.5 m	3.0 m						
9903535	3.5 m	3.5 m						

Measure 3 : This measurement is generally 0.5 m less than measurement 2

Maximum capacity of the rack	
KWIKtrak 180 mm	
kg	lb
270	600



√ = Included / o = Optional

# Part to be ordered separately	Description		
a	700.10650	Trolley - KWIKtrak (for portable ceiling lifts)	o
b	700.15510	Charger interface (for fixed ceiling lifts)	o
c1	700.08310	Reacher - 61 cm (24 in.)	o
c2	700.08320	Reacher - 91 cm (36 in.)	o



Portable non-permanent track system


Semi-permanent rack

Our KWIKtrak semi-permanent rack is a robust free standing system with the capacity to handle everything including your bariatric lifting tasks

Suitable for use with fixed ceiling lifts or one of our portable lifts. Our semi-permanent rack is the ideal solution for temporary yet portable heavy duty lifting.

The semi permanent rack is available in a 2 post or 4 post (X-Y) configuration.

Depending on the choice of track and configuration, it can offer a lifting capacity of up to 1000 lb (455 kg).



# Product	Description
700.11900	Semi-permanent rack KWIKtrak

√ = Included / o = Optional

# Part	Description	
1	700.10650 Trolley- KWIKtrak (for portable ceiling lift)	o
2	200.12810 KWIKtrak 90 mm - Straight 3.0 m	o
3a	200.12872 KWIKtrak 140 mm - Straight 3.5 m	o
3b	200.12870 KWIKtrak 140 mm - Straight 4.5 m	o
4a	200.12882 KWIKtrak 180 mm - Straight 3.5 m	o
4b	200.12884 KWIKtrak 180 mm - Straight 4.5 m	o
4c	200.12880 KWIKtrak 180 mm - Straight 6.0 m	o
5a	700.08310 Reacher - 61 cm (24 in.)	o
5b	700.08320 Reacher - 91 cm (36 in.)	o

Maximum KWIKtrak length vs capacity and track size			
Lift capacity	KWIKtrak		
	90 mm	140 mm	180 mm
200 kg (440 lb)	3 m	4.5 m	6 m
200 kg (440 lb)		3.5 m	4.5 m
455 kg (1000 lb)			3.5 m

