

Everyday Electric Lift

- Heavy-gauge steel construction
- Six-point spreader bar with 360° rotation
- Uses Lumex® six-point, four-point slings or 2-point slings
- Manually adjustable base
- 24V DC motor allows quiet smooth operation
- Ergonomically designed hand pendant with pendant hook
- Easy access emergency stop button
- Maximum weight capacity 400 lb
- Optional digital scale available
- Designed to meet the requirements of HCPCS code:E0635
- 1 year limited warranty

LF1040

Blue

LF1041

Black/Silver Hammertone

Specifications

Type.....	Electric	Overall Width.....	25.6"
Max. Weight Capacity.....	400 lb	Max. Base Width.....	39.4"
Floor Clearance.....	5.5"	Min. Base Width.....	25.6"
Lifting Range Min. Height.....	27.4"	Weight.....	115.0 lb
Max. Height.....	77.7"	Motor.....	24v DC
Overall Length.....	45.3"	Warranty Frame.....	1 year
Overall Height.....	55.1"	Battery.....	1 year



GRAHAM-FIELD

www.grahamfield.com