



Joerns Healthcare is committed to providing a complete line of top quality lifting and repositioning equipment to reduce workplace injury risks and overall facility insurance costs. Our Hoyerpro+ lifts have been designed to work together to *Create a Safer Environment (C.A.S.E)* for your staff and patients, while addressing your budget requirements.

The Joerns Hoyerpro+ 600 Ceiling Lift raises ceiling lift functionality and usability to a new level. The compact design and modern styling allows it to blend in with any décor. Small and sleek, this lift can safely handle up to 600lbs. Designed specifically for the institutional market, the HoyerPro+ 600 Ceiling Lift features a high output lifting motor and an all-metal drive train for superior safety and durability. Our 6-point spreader bar means maximum comfort and safety for your patients. It also allows the lift to be used for effective in-bed repositioning. The Hoyerpro+ 600 Ceiling Lift is available in both manual and power traverse models, and is packed with an extensive array of safety and ease-of-use features.

The Hoyerpro+ 600 Ceiling Lift is guaranteed to reduce your workplace injuries and increase patient safety. In addition, it is the only ceiling lift that can be seamlessly integrated into your environment for the best in convenience, appearance and patient privacy.

Joerns® Hoyerpro+ 600 Ceiling Lift Features

- Industry exclusive trolley design accommodates curtain track solution.
- Up to 600lb lifting capacity. Choose manual or power traverse models.
- Compact design with high output motor design.
- Integrated safety features include slack strap cutout switch, angle detector, anti-freewheel, limit switches, overload protection, emergency mechanical lowering and integrated cycle counter.
- Extra wide and extra long lifting strap for safety and convenience.
- Super six point spreader bar for superior safety and comfort.
- Emergency on/off.
- Two or Four function hand controller electric operation for optimal reliability.
- Smart sleep mode function to conserve battery power.
- Digital control pad with visual and audio low battery indicators.
- All metal lifting and traverse drive train for safety, reliability and durability.

Options

Digital Weigh Scale



The optional digital weigh scale allows for accurate patient weighing in lbs. and kgs.

Curtain Track Pass Through



A patent pending trolley and track design allows your existing privacy curtain tracks to work seamlessly with ceiling lift tracks.

You must use Hoyer® slings with Hoyer lifts.
Specifications may vary due to manufacturing tolerances.

Specifications

Power Source:	2x12V SLA Batteries, 4.5 AH
Safe Working Load	
Model 600:	600 lbs (272 kgs)
Width:	8" (20.3 cm)
Length:	16" (30.64 cm)
Spreader Bar Length:	23" (58.42 cm)
Lifting Strap Width:	2" (36.6 cm)
Lifting Strap Length:	8' (243 cm)
Charger:	120VAC/24VDC CSA/UL Approved
Lifting Motor:	24VDC
Lifting Motor Drive System:	Direct Drive Gear
Traverse Motor (T7003):	24 VDC
Traverse Motor Drive System:	Chain/Sprocket
Overall Weight:	30 lbs (13.6 kg)
Charging:	End of Track
Emergency Stop:	Standard
Visual Battery Indicator:	Standard
Audio Low-Battery Indicator:	Standard
Emergency Lowering:	Standard
Slack Strap Detector:	Standard
Angle Detector Cutout Switch:	Standard
Anti Free Wheel Safety System:	Standard
Up Limit Switch:	Standard
Down Limit Switch:	Standard
Remote hand Controller:	Standard
Warranty:	1 Year
Digital Weigh Scale:	Optional
Approved:	CSA-Z10535-03



You must use Hoyer® slings with Hoyer lifts.
Specifications may vary due to manufacturing tolerances.

