

WE'RE THE COMPANY THAT BRINGS YOU FGSH 450

HYDRAULIC FLOOR LIFT



Open



Closed



Our new hydraulic FGSH 450 lift was created to make transfer situations safer and more affordable. This manual lift is the perfect combination of functionality and economic value. It is versatile enough to use for most patient transfers. The features include; a padded swivel bar and push handle, 450 lbs weight capacity, durable steel construction with an attractive painted finish, and a 360 degree carry bar rotation without side-to-side sway. Lightweight construction and easy disassembly allow for quick transport and setup. Most importantly, the FGSH 450 offers safety, comfort, and stability to any Safe Patient Handling program and meets the challenges faced in a homecare setting.

SPECIFICATIONS

- Base width: 22" - 42.25"
- Base height: 6.5"
- Base length: 46.75"
- Boom height at sling attachment: 17" - 61"
- Boom overall height: 73"
- Weight limit: 450 lbs
- Caster size: 5"

PRODUCT HIGHLIGHTS

- Six-point swivel bar attachment easily adapts to all sling styles and positions
- Large casters significantly improve navigational control and contribute to caregiver safety
- 17" - 61" range allows for patient to be picked up from supine position on the floor
- Pump handle can rotate from side to side for convenience of attendant
- Carry bar is padded for extra patient safety
- Completely enclosed mechanism controls adjustable base and allows for easy cleaning and disinfection
- Base offers obstruction-free operation and convenient locking mechanism to hold legs in position

PRISM MEDICAL

Prism Medical

45 Progress Parkway, Maryland Heights, MO USA 63043 | T 314.692.9135 | 1.866.891.6502

www.prismmedicalinc.com