



# XO-404

## Horizon Standing Wheelchair Manual Propel

Utilizing the lightest and strongest aluminum in the world, as well as incorporating ergonomics features into the design, the XO-404 Horizon Standing Wheelchair is the ideal chair for users who are active. The standing wheelchair allows users to independently propel the wheelchair throughout daily activities and manually stand when desired.

With its armrest frame curving to allow for easier transferring and having the utmost appeal in the manual standing wheelchair market, this standing wheelchair has all the ideal features for users to move independently and freely.

# XO-404

Horizon Standing Wheelchair Manual Propel

## PRODUCT FEATURES



- 57 lbs total weight with wheels
- 6061 T-6 Aircraft Grade Aluminium Frame
- Pure white frame color
- Transparent lightweight composite side panel
- High strength struts for manual standing
- Rear strut for stability and increased nimbleness
- Easy to use standing handle allows for manual stand in as little as 10 seconds
- Quick release detachable rear wheels
- Lightweight frame design
- Breathable nylon sling/upholstery seat and back
- 24" x 1" high profile pneumatic tires
- 6" x 1 1/4" poly foam front flat free casters
- Padded flip-back armrest
- Flip up footplates
- Chest strap, pelvic seatbelt, and knee straps are standard
- Scissor style wheel
- Lifetime warranty on frames

## SPECIFICATIONS

HCPCS: N/A

Overall Height: 32 inch

Overall Width: 29 inch

Overall Length: 42 inch

Armrest Height: 8 inch

Seat Width: 18 inch

Seat Depth: 18/19/20 inch

Seat Height: 21 inch

Back Height: 11 inch

Weight Capacity: 250 lbs

