

LT-770Q Red Streak

Lightweight Compact Wheelchair



- Detachable/Padded desk length armrest with alloy side guard
- Swing-away/Length adjustable footrest with heel loop footplate
- Push to lock wheel brakes
- High-grade/Padded flame retardant nylon upholstery
- 24" x 1" Quick release Mag wheels with polyurethane tires
- Tool free removable 8"x1" front casters
- Pouch on backrest allowing storage space for items
- Frame color: Red Streak
- Weight capacity: 300 lbs
- Limited Lifetime Warranty

SPECIFICATIONS (subject to change without notice) :

	seat width	seat depth	armrest height	seat height	back height	overall height	overall open width	folded width	overall length	weight w/o riggings	weight capacity
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
LT-770Q	18	16	9	19	17	36	26	12	43	37	300

OPTIONAL ACCESSORIES:

- Elevating Legrests - Anti-Tippers - Brake Extension - Seat Belt - Leg Straps - Headrest - Carrying Pouches – Seat Cushions
- Back Cushions



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335