

# LT-950

*aluminum high strength ultra lightweight, super compact wheelchairs*



- Ultra light weight aluminum frame
- Medicare code: K0004
- Folding backrests
- Fixed armrests & footrests with velcro mount leg strap
- Flat free front casters
- 22" flat free rear wheels with polyurethane tires
- Breathable nylon upholstery
- Light weight only: 25 lbs w/footrests
- Weight capacity: 250 lbs
- Frame color: Black

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	Weight w/riggings (lb.)	weight capacity (lb.)
LT-950	18	15	19 1/2	19	37	26	25	250

Model No.	Description
LT-950	18" seat width, aluminum, high strength ultra lightweight

Common Options: • Seat Belt • Desk Tray

\* Price and specifications are subject to change without notice.



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418  
 Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335