

S-Ergo 106


Lightweight Ergonomic Manual

Wheelchair W/ User Angle Adjustable Back



Features for S-Ergonomic Series

- Ergonomic ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slippage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider/Concaved arm pad, providing comfort & support



- Medicare code K0004
- S-Style Ergonomic Seat
- High strength, weighs only 27 lbs. (w/ footrests)
- Silver ¼” Anti-Bacterial cushion, able to be removed to wash
- Has 3 easy angle adjustments: “Working” Angle (95°), “Moving” Angle (100°) & “Resting” Angle (105°)
- “Rigidifying” push-bar & folds to strengthen back support, comes standard on this model!
- Folding backrest/folding seat for transporting vehicle or travel
- Extending push to lock manual breaks
- Fixed armrests w/wider concaved armpads
- Fixed footrests
- “Tube-in Center” foot-plate, assures side leg support
- Pouch for carrying small items attached to upholstery behind backrest.
- 24” Rear polyurethane, high tread, flat free wheels
- 7x1” Front Polyurethane casters
- 16x17, 18x17 Seat widths available
- Anti-Tippers standard on this model!
- Frame Color: Pearl Silver Only
- Upholstery: Black breathable mesh bottom/ AEGIS top

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
S-106F16	16	17	8	19	17	36	23	12	39	27	250
S-106F18	18	17	8	19	17	36	25	12	39	27	250

**Available Options ~ Seat belt, Memory Foam Cushion, Adjustable UV Resistant Canopy, Rear Wheel Reflector/Spoke Guard, Magnesium Wheels, Vulcan Grip Wheels, Foldable/ Adjustable Headrest, Frog Leg Suspension*



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry. CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335