



## T-900W Series

### *Deluxe Steel Bariatric Transport Chair*

- ❖ Padded detachable desk length armrests
- ❖ 19" seat with carry loops for transporting convenience
- ❖ Swing-away footrests w/ heel loop
- ❖ 8" front caster, 12 1/2" rear solid polyurethane tire
- ❖ Manual lock brakes
- ❖ High-grade flame retardant, breathable nylon upholstery
- ❖ Dual Cross Bars
- ❖ 450 Weight Capacity
- ❖ Red Frame Color

Model No	Seat Width (in.)	Seat Depth (in.)	Seat Height (in.)	Back Height (in.)	Overall Height (in.)	Overall Open width (in.)	Overall Weight (lb.)	Weight Capacity (lb.)
T-920W	20	18	21	16	36	25	40	450
T-922W	22	18	21	16	36	27	44	450

\*\* Price & specifications are subject to change without notice.



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue. City of Industry. CA 91748-1418  
 Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335