

GO CHAIR®

Get Out There!



Dual under-seat compartments for convenient storage

Shown in:
Ruby Red



The Go Chair® is re-engineered from the ground up, offering a bold new style available in an array of contemporary colors. Enhanced performance and comfort, along with feather-touch disassembly, allows you to enjoy lightweight travel and independence on the go. With an increased weight capacity of 300 lbs. and a sleek, bold look, the Go Chair® makes travel easy.



One-hand feather-touch disassembly

Sapphire Blue Rose Quartz Citrine Yellow Amber Orange Pearl White



Variety of color options

SPECIFICATIONS

GO CHAIR®

Drive Wheels:	9" (solid)
Front Casters:	6" (solid)
Rear Anti-tips:	3" (solid)
Maximum Speed ¹ :	Up to 3.7 mph
Ground Clearance ² :	1.6" at anti-tip bracket 1.9" at frame
Turning Radius ² :	25.4" with foot platform
Length:	33" with front rigging
Width ⁶ :	22.6"
Seat Type:	Foldable, compact, lightweight
Seat to Ground Height:	22"-23"
Batteries ^{3,5} :	(2) 18 Ah, 12 V, deep-cycle
Range ^{1,4} :	Up to 8.74 miles
Charger:	Off-board, 2 amp
Electronics:	40 amp, PG GC3 (non-programmable)
Weight Capacity:	300 lbs.
Base Weight:	67 lbs. (front & rear)
Standard Seat Weight:	30 lbs.
Battery Pack Weight ⁵ :	31 lbs.
Suspension:	Limited
Heaviest Piece:	36 lbs. (rear)

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition tire type and condition. These specifications can be subject to a variance of +10%, -5%.
² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
³ AGM or Gel-Cell type required. *See "Batteries and Charging."
⁴ Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
⁵ Battery weight may vary based on manufacturer.
⁶ Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA and ISO 7176 series standards. This product is not coded by DME PDAC or does not meet criteria: K0899. All specifications subject to change without notice.

LIFT FINDER

Mobility Made Easy. Pride offers a variety of lifts to transport your mobility product with ease. To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at pridemobility.com/lift-finder/

4-wheel

FEATURES

- Choice of 6 colors
- 18" x 17", 5-point swivel seat with latch release
- Improved agility - longer wheel base; wider stance; compact turning radius 25.4"
- Dual swing-away, under seat storage bins
- Height and width adjustable armrests
- Larger footplate
- One hand feather-touch disassembly
- 60" lap belt

OPTIONS

- Cup holder
- Lap belt- 50", 70"
- Safety flag
- Oxygen holder
- Dual saddle bag
- Weather cover
- Additional battery pack
- Rear mount bracket cane/crutch holder
 - Single (3")
 - Double (3")
- Armrest mount cane clip

*Military Patches Available



Pride
Mobility Products Corp.
— LIVE YOUR BEST —

GO GO
TRAVEL MOBILITY

US: 800-800-8586 • CAN: 888-570-1113

pridemobility.com