

Jazzy®

600 ES

The Jazzy® 600 ES features patented Mid-Wheel 6® technology and Active-Trac® ATX Suspension for an extraordinary level of performance, maneuverability and the tightest turning radius in the Jazzy line. Mid-Wheel 6 technology incorporates ATX Suspension for extra stability with independent front casters that raise and lower as needed to compensate instantly for changes in terrain and ease of transitions. The Jazzy 600 ES offers true all-around performance indoors and out, transitions curbs smoothly and is rated for climbing and aggressive performance, making it an excellent power chair choice.

*Jazzy® 600 ES
shown in Cardinal Red
with high-back seat*



Pride®
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com



SPECIFICATIONS

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	2.75"
Turning Radius	20.5"
Overall Length	35.5"
Overall Width	24.5"
Drive Wheels	14"
Front Wheels	6" solid OMNI-casters
Rear Wheels	6" solid OMNI-casters
Suspension	Active-Trac® (ATX)
Drivetrain	Two-motor, mid-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	50A, Dynamic Shark controller
Battery Charger	Off-board, 5 amp
Per-Charge Range	Up to 16.65 miles
Battery Requirements	(2) 12 volt, deep cycle, NF-22
Battery Weight	37 lbs. each
Base Weight	123.6 lbs.
Standard Seat Weight	51.5 lbs. (high-back)
Standard Seat	18" x 18-20" and 20" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors

COLOR CHOICES



Cardinal
Red

OPTIONAL SEATING

- Synergy® Seat
- Solid seat pan

MODELS

Jazzy 600 ES 2S-SS	Code: K0822
Jazzy 600 ES 2S-C	Code: K0823

FEATURES

- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Heel loops
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Military patches
- Swing-away joystick
- XLR USB charger

RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Olympian
- Commander

LIFT FINDER

Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com

Pride®
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.