


AL560 Automatic Universal Power Chair Lift



FEATURES

- Carries many popular power chairs
- Padded hold-down arm conveniently secures power chair
- No power chair modifications required
- Platform adjusts to provide a snug, custom fit
- Manual crank backup & license plate mount included
- Swing-away option available
- Lift folds automatically when unoccupied

SPECIFICATIONS

Lifting Capacity:	350 lbs (160 kgs)
Platform Size (Useable Space):	28.5" x 38.5", adjustable
Installed Weight:	105 lbs (48 kgs)
Hitch Height (Range):	5" - 19"
Hitch Class:	II or III
Warranty:	3-year transferable

DESCRIPTION

Harmar's AL560 is designed to meet the needs of consumers who need a quick, automatic securing option without making modifications to their power chair.

The adjustable platform and vertically adjustable, padded hold down arm will safely and automatically secure virtually any power chair up to 350 lbs. When empty, the platform automatically folds up and out of the way. One switch operation and drive on/off platform make this fully automatic, American Made lift even easier to use!