

PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

# GO GO ELITE TRAVELLER®

3 & 4-WHEEL SCOOTER

THE  
QUICKEST WAY  
TO TRAVEL



Feather-touch disassembly

The compact design of the **Go-Go Elite Traveller®** allows you to easily maneuver in tight spaces while providing stable outdoor performance. With a 300 lb. weight capacity, a charger port in the tiller, one-hand feather-touch disassembly and maximum speeds up to 4 mph make it easy to go where you want to go with the Go-Go Elite Traveller.



Go-Go Elite Traveller®  
3-wheel



Go-Go Elite Traveller®  
4-wheel

[pridemobility.com](http://pridemobility.com)



# SPECIFICATIONS

## GO-GO ELITE TRAVELLER® 3W: SC40E 4W: SC44E

Weight capacity:	300 lbs.
Maximum speed <sup>1</sup> :	Up to 4 mph
Ground clearance <sup>2</sup> :	2"
Turning radius:	3W: 33" 4W: 44"
Overall length <sup>2</sup> :	3W: 37" 4W: 39.5"
Overall width <sup>2</sup> :	19.5"
Front tires:	3W: 2" x 8" solid 4W: 2" x 7" solid
Rear tires:	2.5" x 8" solid
Range per charge <sup>1,4</sup> :	3W: Up to 7.2 miles (12 AH) 4W: Up to 6 miles (12 AH) 3W: Up to 10.4 miles (18 AH) 4W: Up to 9.7 miles (18 AH)
Total weight without batteries <sup>5</sup> :	3W: 76 lbs. 4W: 80 lbs.
Heaviest piece when disassembled:	3W: 29.5 lbs. (front section) 4W: 34 lbs. (front section)
Standard seating:	Type: Foldable, compact, vinyl-covered molded plastic Weight: 23 lbs. Material: Black Vinyl Dimensions: Width: 17", Depth: 17"
Motor output:	24-volt DC motor, sealed mini-transaxle
Brakes:	Regenerative and electromechanical
Battery requirements <sup>3,5</sup> :	(2) 12V Deep-cycle Size: (2) 12 AH, 9 lbs. each or (2) 18 AH, 14 lbs. each (optional) Weight: 12AH: 20 lbs. per battery pack, 18 AH: 29.5 lbs. per battery pack
Battery charger:	Off-board, 2 amp
Color:	2 sets of exciting colors with every scooter!



Red



Blue

- <sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- <sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- <sup>3</sup> AGM or gel cell type required.
- <sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- <sup>5</sup> Battery weight may vary based on manufacturer.

**NOTE: All specifications subject to change without notice.**

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

## FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand easily disassembles into five super lightweight pieces for convenient transport and storage
- Pride's exclusive black, non-scuffing tires
- Includes sets of easily changeable colored panels - red and blue
- Front frame-mounted seat post offers maximum stability
- Microprocessor-based controller offers optimal power management and added safety features
- Standard front tiller-mounted basket
- Removable footboard basket for additional storage (4-wheel version only)
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board
- Charger port in tiller

## ACCESSORIES

- Forearm crutch holder
- Rearview mirror
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Quad cane holder
- Wishbone crutch holder
- Double crutch holder
- Weather cover
- Saddlebag
- Cell phone holder
- Cane clip attachment
- Lap belt
- Under seat storage (assembly required)
- Additional 12AH battery box
- Additional 18AH battery box



## LIFT FINDER

**Mobility Made Easy.** Pride offers a variety of lifts to transport your mobility product with ease. To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at [pridemobility.com/lift-finder/](http://pridemobility.com/lift-finder/)