

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

GO-GO® ULTRA X

3 & 4-WHEEL SCOOTER

ECONOMY CLASS PRICE, FIRST CLASS FEATURES

The Go-Go® Ultra X is loaded with features such as one-hand disassembly and a convenient drop-in battery box for worry-free travel. With a 300 lb. weight capacity, a maximum speed up to 4 mph, a per-charge range up to 10.4 miles on the 3-wheel and 13 miles on the 4-wheel, Go-Go® Ultra X is an exceptional value.



Go-Go® Ultra X 3-wheel



Go-Go® Ultra X 4-wheel




Feather-touch disassembly

pridemobility.com



SPECIFICATIONS

GO-GO® ULTRA X 3W: S40X 4W: S44X

Weight Capacity:	300 lbs.
Maximum Speed ¹ :	Up to 4 mph
Ground Clearance ² :	2.25" at mid-frame, 1.25" at motor
Turning Radius ² :	3W: 33" 4W: 54.43"
Overall Length ² :	3W: 37" 4W: 40.16"
Overall Width ² :	19.5"
Front Tires:	3W: 2" x 8" solid 4W: 2" x 7" solid
Rear Tires:	2.5" x 8" solid
Range Per Charge ^{1,4} :	3W: 6.9 miles (12 AH) 4W: 8.2 miles (12 AH) 3W: 10.4 miles (18 AH) 4W: 13.1 miles (18 AH)
Total weight without batteries ⁴ :	3W: 76 lbs. 4W: 80 lbs.
Heaviest Piece When Disassembled:	3W: 29.5 lbs. (front section) 4W: 34.4 lbs. (front section)
Standard Seating:	Type: Foldable, compact, vinyl-covered molded plastic Weight: 23 lbs. Material: black vinyl Dimensions: Width: 17", Depth: 17"
Motor Output:	24-volt DC motor, sealed mini-transaxle
Brakes:	Regenerative and electromechanical
Battery Requirements ^{2,3} :	(2) 12V deep-cycle Size: (2) 12 AH (standard) or 18 AH (optional) Weight: 9 lbs. each, 20 lbs. per battery pack, 14 lbs. each, 29.5 lbs. per battery pack
Battery Charger:	Off-board, 2 amp
Color:	 Red

FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Black, non-scuffing tires exclusively from Pride®
- Front frame-mounted seat post offers maximum stability
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board

ACCESSORIES

- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Quad cane holder
- Double crutch holder
- Weather cover
- Saddlebag
- Underseat storage
- Cane clip attachment
- Optional battery pack
- Under seat storage (assembly required)



LIFT FINDER

Mobility Made Easy. Pride offers a variety of lifts to transport your mobility product with ease. To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at pridemobility.com/lift-finder/

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³ AGM or gel cell type required.

⁴ This range was obtained with the 4-wheel model.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.