

R82[®] 



toucan[™]



A standing frame to increase quality of life
happiness

The Toucan standing frame appears with soft curves and smart colours. It is designed to support the user in the upright position.


The distance between the left and the right part of the frame allows almost any wheelchair to go in between. A good feature for transfer from wheelchairs ect. into the Toucan.

The Toucan is available in 3 sizes and suitable for children from the age of aprox. 1-10 years. Not all children are alike which is why we recommend to use the measuring chart to find the correct size for your child.


toucan™ by **R82**®

Measurements


	Size 1 cm (inch)	Size 2 cm (inch)	Size 3 cm (inch)
Length (L)	75 (29 1/4")	88 (34 1/4")	88 (34 1/4")
Width (W)	56 (21 3/4")	63 (24 1/2")	63 (24 1/2")
Total height (H1)	67-87 (26 1/4-34")	72-101 (28-39 1/2")	97-123 (37 3/4-48")
Height from foot/chest (H2)	34-70 (13 1/4-27 1/4")	54-90 (21-35")	74-110 (29-43")
From floor to foot plate	6 (2 1/2")	6 (2 1/2")	6 (2 1/2")
Angle of the centre bar (A1)	0 to +30°	0 to +30°	0 to +30°
Angle of the foot plate (A2)	+15° to -10°	+15° to -10°	+15° to -10°
Weight	7,5 kg (16,5 lb)	8,5 kg (19 lb)	9 kg (20 lb)
Max. load/user weight	30 kg (66 lb)	40 kg (88 lb)	50 kg (110 lb)



Size 1



Size 2+3



The Toucan standing frame appears with soft curves and smart colours. It is designed to support the user in the upright position.

The distance between the left and the right part of the frame allows almost any wheelchair to go in between. A good feature for transfer from wheelchairs etc. into the Toucan.

The Toucan is available in 5 sizes and suitable for children/teenagers from the age of approx. 1-18 years. We recommend to use the measuring chart below to find the correct size.




Knee support, hinged



Tray with bowl & cover plate



Hinged sandals



Measurements

	Size 1 cm (inch)	Size 2 cm (inch)	Size 3 cm (inch)	Size 4 cm (inch)	Size 5 cm (inch)
Length (L)	75 (29 ¹ / ₄ "	88 (34 ¹ / ₄ "	88 (34 ¹ / ₄ "	111 (43 ¹ / ₄ "	111 (43 ¹ / ₄ "
Width (W)	56 (21 ³ / ₄ "	63 (24 ¹ / ₂ "	63 (24 ¹ / ₂ "	70 (27 ¹ / ₄ "	70 (27 ¹ / ₄ "
Total height (H1)	67-87 (26 ¹ / ₄ -34"	72-101 (28-39 ¹ / ₂ "	97-123 (37 ³ / ₄ -48"	104-134 (40 ¹ / ₂ -52 ¹ / ₄ "	124-154 (48 ¹ / ₄ -60"
Height from foot/chest (H2)	34-70 (13 ¹ / ₄ -27 ¹ / ₄ "	54-90 (21-35"	74-110 (29-43"	91-121 (35 ¹ / ₂ -47 ¹ / ₄ "	111-141 (43 ¹ / ₄ -55"
From floor to foot plate	6 (2 ¹ / ₂ "	6 (2 ¹ / ₂ "	6 (2 ¹ / ₂ "	6 (2 ¹ / ₂ "	6 (2 ¹ / ₂ "
Angle of the centre bar (A1)	0 to +30°	0 to +30°	0 to +30°	0 to +20°	0 to +20°
Angle of the foot plate (A2)	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°
Weight	7,5 kg (16,5 lb)	8,5 kg (19 lb)	9 kg (20 lb)	10 ¹ / ₂ kg (23 lb)	10 ³ / ₄ kg (24 lb)
Max. load/user weight	30 kg (66 lb)	40 kg (88 lb)	50 kg (110 lb)	60 kg (132 lb)	70 kg (154 lb)

R82 A/S

Parallevej 3
8751 Gedved
Denmark
Tel +45 79 68 58 88
Fax +45 75 66 51 92
E-mail: r82@r82.com
Web: www.r82.com

R82 UK Ltd.

Unit D4A
Coombswood Business Park East
Coombswood Way, Halesowen
West Midlands B62 8BH
Great Britain
Tel 0121 561 2222
Fax 0121 559 5437
E-mail: info@r82.com
Web: www.r82.com

Distributor: