

SPECIFICATIONS:

Pump:

- Power: 120 VAC
- Input Frequency: 60 Hz
- Pressure Range: 20 mmHg to 60 mmHg (inner mattress pressure)
- Auto firm (60 mmHg, Auto return in 30 min.)
- Seat function (Pressure preset + 10 mmHg)
- Static Selectable
- Static auto return (in 2 hours)
- Comfort therapy
- Function lock (Auto lock in 30 min. after on)
- Alarm mute
- Low pressure alarm
- System (alternating) passive alarm
- Power failure alarm
- Low pressure alarm delay
- Memory function
- Air Flow: 12-LPM (Litres Per Minute)
- Functions: Static and dynamic features
- Power Cord Length: 15 Ft.
- Pump Dimensions: 14.2" (W) x 9.8" (H) x 4.5" (L)
- UL: 60601-1, 1st Edition-2003
- Maximum Current: 0.5 Amps
- Warranty: 18-Months Limited

Mattress:

- Topper Material: Zippered and Quilted Poly/PU Urethane Fabric (4-Way Stretch) is; water resistant, vapor permeable, machine washable and low shear
- Base Material: Zippered Vinyl
- 4-Air cells at head section remain static to enhance pillow function
- Dimensions: 36" (W) x 80" (L) x 8" (H)
- Maximum Support Capacity: 450 lbs.
- Warranty: 18-Months Limited