

## **SPECIFICATIONS:**

### **Pump:**

- Power: 120 VAC
- Input Frequency: 60 Hz
- Pressure Range: 30 mmHg to 80 mmHg (inner mattress pressure)
- Air Flow: 8-LPM (Litres Per Minute)
- Selectable Cycle Times: 10, 15, 20 and 25-Minutes
- Functions: Static and dynamic features
- Power Cord Length: 15-Ft.
- Pump Dimensions: 13.1" (W) x 9.7" (H) x 4.7" (L)
- UL: 60601-1, 1st Edition-2003
- Maximum Current: 0.5 Amps
- Warranty: 18-Months Limited

### **Mattress Dimension:**

- Topper Material: Zippered and Quilted Urethane Coated Nylon Fabric (4-Way Stretch) is; water resistant, vapor permeable, machine washable and low sheer
- Base Material: Zippered Vinyl
- 4-Air cells at head section remain static to enhance pillow function
- Dimensions: 60" (W) x 80" (L) x 8" (H)
- Maximum Support Capacity: 1000 lbs.
- Warranty: 18-Months Limited