

COMPARE TO:

JAY® BASIC BACK



14889

### General Use Back Cushion with Lumbar Support

- 1P** 14906 16" (W) x 17" (D), 1/ea
- 1P** 14889 18" (W) x 17" (D), 1/ea
- 1P** 14920 20" (W) x 17" (D), 1/ea

HCPCS code: E2611

**FEATURES:**

- Improves seating posture while providing back support and comfort.
- **Durable Masonite** board secured on back of foam core prevents hammocking in sling-back wheelchairs.
- Conforms to lumbar curve to support lower back.
- Sealed liner protects foam core from moisture.
- Has a non-slip vinyl back to prevent the cushion from sliding.
- 300 lb. weight capacity.
- 18 Month Limited Warranty on cover and cushion.

### Premium Guard Gel Overlay

- 1P** 14893 Gel Overlay, 76" (L) x 34" (W) x 3.5" (H), 1/ea
- 1P** 14901 Bariatric Gel Overlay, 76" (L) x 42" (W) x 3.5" (H), 1/ea

HCPCS code: E0185



14893

**PREMIUM GUARD GEL OVERLAY WITH SOFT NYLON COVER!**

**300 LB**  
WEIGHT CAPACITY

**FEATURES:**

- New and improved gel overlay has an upgraded **300 lb. weight capacity** and optimizes pressure redistribution.
- **New upgrade: 15 gel bladders cover entire patient contact area.**
- Full length "head-to-toe" gel bladders minimizes "bottoming out" and "pitting".
- New convoluted foam topper vastly improves patient comfort and air circulation.
- 4 elastic corner straps secures overlay to mattress.
- Urethane coated nylon cover is water resistant and vapor permeable.
- Vinyl base is stain resistant and waterproof.
- Bladders contain an aqueous non-toxic gel and polystyrene beads.
- Designated head and foot sections.
- 300 lb. weight capacity (14893), 500 lb. weight capacity (14901).
- 18 Month Limited Warranty.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

**1P = 1-Pack Available**

**SPECIFICATIONS**

Item #	Width	Depth	Height	Product Weight	Carton Shipping Weight	Weight Limit
14906	16"	17"	2.5"	2 lbs.	3.0 lbs.	300 lbs.
14889	18"	17"	2.5"	2.7 lbs.	3.0 lbs.	300 lbs.
14920	20"	17"	2.5"	3.7 lbs.	4.25 lbs.	300 lbs.
14893	76"	34"	3.5"	47.8 lbs.	50.8 lbs.	300 lbs.
14901	76"	42"	3.5"	50.5 lbs.	53.5 lbs.	500 lbs.