



Two carry handles on side make lifting and placement easier!

NEW!

Non-slip bottom and 4 corner straps keep the mattress securely in place! (underside shown)

Available 1 Pack Quick 'N Easy Comfort Mattress
 15076 80" (L) x 36" (W) x 6" (H), 1/ea

HCPCS Code: E0272

FEATURES:

- Premium quality homecare mattress offering the utmost in comfort.
- **Patent pending super compression technology packages mattress in a box with a footprint of just 14" x 14" x 37.5".**
- Decompresses quickly.
- Mattress core made from select grade of polyurethane. Offers comfort, feel and advantages of latex rubber, without latex allergen.
- Inner core encased in fire retardant barrier, meeting regulation 16 CFR 1633.
- Non-slip bottom and 4 corner straps keep the mattress securely in place.
- Two carry handles on side make lifting and placement easier.
- Durable nylon cover is waterproof for easy use and care.
- 275 lb. weight capacity.
- Limited 1 Year Warranty.

SPECIFICATIONS

Item #	Length	Width	Height	Product Weight	Carton Shipping Weight	Weight Capacity
15076	80"	36"	6"	20 lbs.	22 lbs.	275 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.