



Gravity 7 Standard Long Term Care Pressure Redistribution Mattress

The Gravity 7 is an incredible value for a full, 3 layered pressure redistribution foam mattress with a true memory foam heel section. Combined with an intelligently engineered die-cut top layer, and full length articulation cuts on the bottom layer, this product provides more features and benefits, longer warranties and better value than the competition.



3 LAYER DESIGN

**MULTI-LAYERED, INTELLIGENT DESIGN
AND A 7 YEAR GUARANTEE ON OUR MOST
ECONOMICAL FOAM MATTRESS OPTION!**



COMPARE TO:

Direct Supply® Panacea
Invacare® Solace Prevention
Mason® 6500 Lite
Span America® GeoMatt Plus
Blue Chip® Relief Care Pro

CUSTOM SIZES AVAILABLE

15770	80" (L) x 36" (W) x 6" (H), 1/ea
15777	80" (L) x 36" (W) x 6" (H) with Raised Side Rails, 1/ea
15784	84" (L) x 36" (W) x 6" (H), 1/ea
15785	84" (L) x 36" (W) x 6" (H) with Raised Side Rails, 1/ea
15876	76" (L) x 36" (W) x 6" (H), 1/ea
15877	76" (L) x 36" (W) x 6" (H) with Raised Side Rails, 1/ea

HCPCS code: E0184

FEATURES:

- Economical option maximizes value and comfort with 3 full layers of foam available with or without raised side rails.
- Top layer uses unique die-cutting and a true memory foam heel section to optimize zoned pressure redistribution.
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress significantly.
- Cover is fire retardant, fluid proof, tear resistant, low shear, and easy to clean.
- **1 Year Warranty on the cover.**
- **Limited 7 Year Warranty.**

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

= 1-Pack Available