

SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11105N-4	18.5"	25"	13.5"	16.5"	22.5"	16.5"	14.1 lbs.	58 lbs.	350 lbs.
11149-4	20.08"	21.5"	13.5"	14.5"	22.24"	16.34"	13 lbs.	57 lbs.	300 lbs.
11149-1	18"	21.5"	13.5"	14.5"	22.24"	16.34"	13 lbs.	15 lbs.	300 lbs.
11149RB/R-1	18"	21.5"	13.5"	14.5"	22.24"	16.34"	13 lbs.	15 lbs.	300 lbs.