

Steel Drop-Arm Commode with Wheels and Padded Arms

- 11101W-2 Assembled, with Wheels, 2/cs
- 11101WL Caster without Lock for 11101W-2, 2/bx
- 11101WLL Caster with Lock for 11101W-2, 2/bx

NEW IMPROVED ARM!



11101W-2



WHEN ARM IS IN FOLDED POSITION, THE ARM FRAME PROVIDES A "HANDLE" SUPPORT FOR TRANSFERS.



FEATURES:

- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.
- Arm provides additional support while transferring.
- Rust resistant, swivel casters (2 with locks).
- White powder coated steel frame is attractive and easy to maintain.
- Padded arm rests provide extra comfort and support.
- Comes with 7.5 qt. commode bucket with carry handle and splash shield.
- Limited Lifetime Warranty.



Biodegradable Sanitary Bag Commode Liner

RTL12085
12 Commode Liners Per Box, 6bx/cs

See Page 15



ENVIRONMENT FRIENDLY PRODUCT

Recycled materials were used to make this product.



NEW REINFORCED TOOL FREE BACK REST DESIGN.



11112-2

Aluminum Commode with Wheels

- 11112-2 Assembled, 2/cs

FEATURES:

- Non-skid, rust resistant, swivel casters (2 with locks).
- Aluminum construction is lightweight and durable.
- Removable, tool free back.
- Angled back for support and comfort.
- Full length plastic arms for added comfort.
- Comes with a 7.5 qt. commode bucket, splash shield and cover. Splash shield acts as pail holder.
- Limited Lifetime Warranty.

NEW IMPROVED CASTERS!



NON-SKID, RUST RESISTANT, SWIVEL CASTERS. TWO WITH LOCKS!

SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11101W-2	18"	24"	14"	16.5"	24"	18"	22.5 lbs.	48 lbs.	300 lbs.
11112-2	18.5"	24"	14"	16.5"	18.75"	18.75"	9 lbs.	21 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.