



ENVIRONMENT FRIENDLY PRODUCT

Recycled materials were used to make this product.

BARIATRIC
650 LB
WEIGHT CAPACITY



1117N-2

Bariatric Folding Commode

Available
1 Pack

- 1117N-2 Folding, (Has replaced 1117-2 and 1116-2), 2/cs
- 1117N-1 Folding, 1/cs
- RTL10389G Replacement Tip, 1" Grey, (Has replaced 10389G-8), 8bx/cs

FEATURES:

- 650 lb. weight capacity.
- Easily opens and folds.
- Folds flat for convenient storage and transportation.
- Durable plastic 16 1/2" deep snap on seat with lid.
- Plastic armrests for added comfort.
- Removable tool free back.
- Comes complete with 12 qt. commode bucket with carry handle, cover and splash shield.
- Limited Lifetime Warranty.

OPENS & FOLDS IN SECONDS WITH SNAP BUTTONS.

1

2

3

See Page 15
Biodegradable Sanitary Bag Commode Liner
RTL12085
12 Commode Liners
Per Box, 6bx/cs

**TOOL FREE ASSEMBLY
IN 3 EASY STEPS!**



TOOL FREE ASSEMBLY

SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11117N-2	24"	31.5"	13.5"	16.5"	21"	15"	17.5 lbs.	41 lbs.	650 lbs.
11117N-1	24"	31.5"	13.5"	16.5"	21"	15"	17.5 lbs.	20.5 lbs.	650 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.