

BARIATRIC

650 LB

WEIGHT CAPACITY



11130-2

ENVIRONMENT FRIENDLY PRODUCT

Recycled materials were used
to make this product.



Deluxe Bariatric Commode

11130-2	Assembled, 2/cs
11129	Replacement Wing Nut, 2/bx
11130B	Replacement Back Rest, 1/ea
11130S	Replacement Seat, 1/ea
11135TIP	Replacement Tip, 1" Grey, 1/ea

FEATURES:

- 650 lb. weight capacity.
- Large, padded backrest provides patient comfort.
- Durable, heavy duty, grey powder coated steel tubing is sturdy and easy to maintain.
- Large, durable snap-on seat.
- Plastic armrests provide comfort and stability.
- Pail easily slides in from the front of the commode.
- Comes complete with 12 qt. commode bucket and cover that removes from behind.
- Limited Lifetime Warranty.



Biodegradable Sanitary Bag Commode Liner

RTL12085
12 Commode Liners
Per Box, 6bx/cs

See Page 15

ENVIRONMENT FRIENDLY PRODUCT

Recycled materials were used
to make this product.



Deluxe Bariatric Drop-Arm Commode

Available
1 Pack

11135-1	Assembled, 1/cs
11135-S	Replacement Seat, 1/ea
11135TIP	Replacement Tip, 1" Grey, 1/ea
11135B	Replacement Back Support, 1/cs

FEATURES:

- 1000 lb. weight capacity.
- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.
- Accommodates individuals up to 1000 lbs.
- Durable, heavy duty, grey powder coated steel tubing is sturdy and easy to maintain.
- Large, durable snap-on seat.
- Padded armrests provide comfort and stability.
- Pail easily slides in from the front of the commode.
- Comes complete with 12 qt. commode bucket and cover that removes from behind.
- Limited Lifetime Warranty.

BARIATRIC

1000 LB

WEIGHT CAPACITY



11135-1

SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11130-2	23.5"	26"	23"	18.5"	22"	15"	20.5 lbs.	46 lbs.	650 lbs.
11135-1	26"	28"	23"	18.5"	23"	19"	22 lbs.	26 lbs.	1000 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.