



**BARIATRIC**  
**850 LB**  
WEIGHT CAPACITY

### Large Bariatric Drop Arm Commode

Available  
**1 Pack**

- 11138-1 Assembled, 1/cs
- 11138B Replacement Back Rest, 1/ea
- 11138SEAT Replacement Seat, 1/ea
- 10374BK-8 Replacement Tip, 1" Black, (Has replaced 11100B), 8bx/cs

**FEATURES:**

- 850 lb. Weight Capacity.
- Extra Large opening in toilet seat, 8 1/2" x 12".
- Durable heavy duty, powder coated steel construction tubing is sturdy and easy to maintain.
- Easy-to-release drop arm mechanism allows for safe lateral transfer to and from commode.
- Extra large, padded, deep seat provides comfort.
- Padded armrests and backrest provide comfort and support.
- Comes complete with extra large commode bucket and cover that removes from behind.
- Limited Lifetime Warranty.



**BARIATRIC**  
**1000 LB**  
WEIGHT CAPACITY

### Bariatric Extra Large Drop Arm Commode

Available  
**1 Pack**

- 11132-1 Assembled, 1/cs
- 11132SEAT Replacement Seat, 1/ea
- 10374BK-8 Replacement Tip, 1" Black, (Has replaced 11100B), 8bx/cs

**FEATURES:**

- 1000 lb. Weight Capacity.
- Easy-to-release drop arm mechanism allows for safe lateral transfer to and from commode.
- Extra large opening in toilet seat, 8 1/2" x 12".
- Durable heavy duty, powder coated steel construction tubing is sturdy and easy to maintain.
- Padded armrests provide comfort and support.
- Comes complete with extra large commode bucket and cover that removes from behind.
- Limited Lifetime Warranty.

**SPECIFICATIONS**

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11138-1	25.5"	27.5"	24"	21"	18.5"	18.5"	37 lbs.	41 lbs.	850 lbs.
11132-1	36.5"	40"	36"	21"	18.5"	NA	49 lbs.	53 lbs.	1000 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.