

SPECIFICATIONS

Item #	Overall Width Open	Overall Width Closed	Overall Height	Seat To Floor Height	Back Height	Width Seat	Depth Seat	Overall Length (with riggings)	Product Weight	Carton Shipping Weight	Weight Limit
RS185001	24.5"	11"	37.5"	22.5"	15"	17"	18"	42"	51 lbs.	63 lbs.	300 lbs.
RS185002	24.5"	11"	37.5"	22.5"	15"	17"	18"	42"	51 lbs.	63 lbs.	300 lbs.
RS185003	23"	9"	37.5"	22.5"	15"	17"	18"	34.75"	48 lbs.	60 lbs.	300 lbs.
RS185005	24.5"	11"	37.5"	22.5"	15"	17"	18"	42"	51 lbs.	63 lbs.	300 lbs.