



12030



LOCKING KNOB



## Open, Padded Raised Toilet Seat with Four Locking Brackets

Available  
**1 Pack**

- 12030 (Has replaced 12031 with plastic seat), 1/cs
- 12030B Replacement Bracket for 12030, 4/bx

### FEATURES:

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- Padded seat provides extra comfort.
- Open seat promotes patient independence.
- Four adjustable locking clamps provides versatility.
- Adjustable angle and seat height.
- Lightweight and portable.
- Easily attaches to toilet bowl and locks in place with two rear locks.
- No tools required for installation.
- Easy to clean.
- Fits most toilets.
- Limited Lifetime Warranty.



12038-3



12040-3

## Raised Toilet Seat

- 12038-3 Bulk Pack, 1/polybag, 3/cs
- 12040-3 Retail, 3 retail boxes/cs, 3/cs

### FEATURES:

- Designed for individuals who have difficulty sitting down, or getting up from the toilet.
- Lightweight and portable.
- Easily attaches to toilet bowl.
- No tools required for installation.
- Easy to clean.
- Fits most toilets.
- Limited Lifetime Warranty.



**NO TOOLS  
REQUIRED FOR  
INSTALLATION!**

**COMPARE TO:**  
**CAREX®** **MEDLINE®**  
 12040 B310-00 MDS80318RH



**ADDS 4" TO HEIGHT  
OF TOILET SEAT!**



### SPECIFICATIONS

Item #	Width (seat)	Depth (seat)	Seat Height	Product Weight	Carton Shipping Weight	Weight Limit
12030	16"	15.5"	5"	6 lbs.	9 lbs.	300 lbs.
12038-3	15"	15"	4"	2.5 lbs.	8.5 lbs.	300 lbs.
12040-3	15"	15"	4"	2.5 lbs.	10.5 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.