



18" WIDTH BETWEEN ARMS!  
12013

FITS ON ELONGATED TOILETS!



**TOOL FREE REMOVABLE ARMS!**



12402

**Available 1 Pack Premium Plastic Elevated, Regular/Elongated Toilet Seat with Lock**

12013 With Arms, 1/cs

**FEATURES:**

- Designed for individuals who have difficulty sitting down, or getting up from the toilet.
- Arms provide added safety and stability.
- Wide opening in front and back for personal hygiene.
- Heavy-duty molded plastic construction provides additional strength and durability.
- Standard locking mechanism ensures safety.
- Lightweight and portable.
- No tools required for installation.
- Easy to clean.
- Fits most elongated toilets.
- Limited Lifetime Warranty.



**TOOL FREE REMOVABLE ARMS!**

**Available 1 Pack Premium Seat Rizer with Removable Arms**

12402 Standard Toilets, (Will become RTL12402 in 2011), 1/cs  
12403 Elongated Toilets (Has replaced 12401), 1/cs

**FEATURES:**

- Designed for individuals who have difficulty sitting down, or getting up from the toilet.
- Allows individual to use their existing toilet seat and lid.
- Tool Free Arms can be removed or added as needed.
- Arm sleeves and receivers are made of metal for more durability.
- Heavy-duty molded plastic construction provides additional strength and durability.
- Lightweight and portable.
- Easily attaches to toilet bowl.
- No tools required for installation.
- Easy to clean.
- Fits most toilets.
- Limited Lifetime Warranty.



**SPECIFICATIONS**

Item #	Width (seat)	Depth (seat)	Seat Height	Width Between Arms	Product Weight	Carton Shipping Weight	Weight Limit
12013	17"	16.5"	4.5"	18"	4.5 lbs.	6.5 lbs.	300 lbs.
12402	14"	17"	3.5"	19"	3.5 lbs.	4.5 lbs.	300 lbs.
12403	14"	19"	3.5"	19"	3.5 lbs.	4.5 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.