




RTL12026



12014



12062

12063



**COMPARE TO:**

**SAVANNAH®**

**LOCKS SECURELY IN PLACE WITH LOCKING KNOB & CLAMP!**

## Locking Elevated Toilet Seat

Available  
**1 Pack**

RTL12026 Without Arms,  
(Has replaced 12026), 1/cs

### FEATURES:

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- Heavy-Duty molded plastic construction provides additional strength and durability.
- Improved locking device is engineered with larger, heavy-duty "worm screw" and locking plate to ensure a safe, secure fit.
- Lightweight and portable.
- Fits most toilets.
- No tools required for installation.
- Easy to clean.
- Limited Lifetime Warranty.

## Premium Plastic Elevated, Regular/Elongated Toilet Seat with Lock

Available  
**1 Pack**

12014 Without Arms, 1/cs

### FEATURES:

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- Wide opening in front and back for personal hygiene.
- Lightweight and portable.
- Standard locking mechanism ensures safety.
- Fits most elongated toilets.
- No tools required for installation.
- Easy to clean.
- Limited Lifetime Warranty.

## Raised Toilet Seat with/without Lid

Available  
**1 Pack**

- 12062 Without Lid, 2", 1/cs
- 12063 With Lid, 2", 1/cs
- RTL12064 Without Lid, 4", (Has replaced 12064), 1/cs
- 12065 With Lid, 4", 1/cs
- 12066 Without Lid, 6", 1/cs
- 12067 With Lid, 6", 1/cs
- 12062-C Replacement Locking Clamps, 2/bx.

### FEATURES:

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- Two hygiene cutouts provide maximum convenience.
- Designed for use with existing toilet seat.
- Lightweight and portable.
- Easily attaches to toilet bowl and locks in place with two rear locks.
- No tools required for installation.
- Easy to clean.
- Fits most toilets.
- Limited Lifetime Warranty.



### SPECIFICATIONS

Item #	Width (seat)	Depth (seat)	Seat Height	Product Weight	Carton Shipping Weight	Weight Limit
RTL12026	17"	16.5"	5"	4.5 lbs.	6 lbs.	300 lbs.
12014	17"	16.5"	4.5"	5 lbs.	7 lbs.	300 lbs.
12062/3	14"	16"	2"	2/3 lbs.	4/5 lbs.	350 lbs.
12064/5	14"	16"	4"	2/3 lbs.	4/5 lbs.	350 lbs.
RTL12064	14"	16"	4"	2/3 lbs.	4/5 lbs.	350 lbs.
12066/7	14"	16"	6"	3.5/4.5 lbs.	6/7.5 lbs.	350 lbs.

\*Important Notice:  
Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.