



**COMES WITH  
ADJUSTABLE WIDTH  
PADDED ARMS!  
(2 POSITION)**



12001-4

**Available 1 Pack Toilet Safety Frame**

- 12001-4 Assembled, 4/cs
- 12001KD-1 Knock Down, 1/cs
- RTL12000 Knock Down, Arms, (Has replaced 12000-4), 1 per Retail Box/4 Boxes/cs
- 12001BR-1 Toilet Safety Frame Bracket, 1/ea

**FEATURES:**

- Anodized aluminum is sturdy and lightweight.
- Armrests provide additional comfort and support.
- Powder coated aluminum bracket easily attaches frame to toilet.
- Arms are height and width adjustable.
- Limited Lifetime Warranty.



RTL2000/12001KD-1



12001BR-1

**FOR USE WITH  
12001-4 & RTL12000  
TOILET SAFETY FRAMES.**



**SPECIFICATIONS**

Item #	Width Between Arms (smallest)	Width Between Arms (largest)	Depth	Height Seat	Product Weight	Carton Shipping Weight	Weight Limit
12001-4	18.5"	22.75"	32"	28"	4 lbs.	18 lbs.	300 lbs.
12001KD-1	18"	25"	32"	28"	4 lbs.	5 lbs.	300 lbs.
RTL12000	19"	23.50"	32"	28"	4 lbs.	20 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

To order call toll free: 877.224.0946

BATH SAFETY

