





### Floor to ceiling pressure fit system

## Easytrack 2 post

The Easytrack is a floor to ceiling pressure fit system requiring no permanent installation. Other versatile options are available such as 3 post, 4 post, bath tub adaptor and wall bracket configurations.



| # Product | Description        |
|-----------|--------------------|
| 9220010   | Easytrack - 2 post |

√ = Included / o = Optional

| # Part                      | Description                                                         |   |
|-----------------------------|---------------------------------------------------------------------|---|
| 1a                          | 700.17250 Extendable track & standard trolley - 1 piece             | √ |
| 1b                          | 700.10500 Standard post - 2 pieces                                  | √ |
| 1c                          | 700.10230 Blue plastic strip - 4 pieces                             | √ |
| 1d                          | 700.10350 Round foot plate - 2 pieces                               | √ |
| 1e                          | 700.10600 Clip-on level - 1 piece                                   | √ |
| 1f                          | 700.10400 Top plate - 2 pieces                                      | √ |
| 1g                          | 700.10720 Round blue cap - 2 pieces                                 | √ |
|                             | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece | o |
| Optional Protective Bag kit |                                                                     |   |
|                             | 002.10450 Easytrack transport bag - 1 piece                         | o |
| Accessories                 |                                                                     |   |
| 2a                          | 700.08310 Reacher - 61 cm (24 in)                                   | o |
| 2b                          | 700.08320 Reacher - 91 cm (36 in)                                   | o |



Floor to ceiling pressure fit system

Easytrack 3 post with pivoting 3<sup>rd</sup> track

The Easytrack is a floor to ceiling pressure fit system requiring no permanent installation.  
Other versatile options are available such as 2 post, 4 post, bath tub adaptor and wall bracket configurations.



| # Product | Description        |
|-----------|--------------------|
| 9330010   | Easytrack - 3 post |

√ = Included / o = Optional

| # Part | Description                                                                      |   |
|--------|----------------------------------------------------------------------------------|---|
| 1a     | 700.10850 Pivoting extendable track 300 cm (118 in) & standard trolley - 1 piece | √ |
| 1b     | 700.10800 Extendable track 300 cm (118 in) & swivel trolley - 1 piece            | √ |
| 1c     | 700.10455 Pivot post - 1 piece                                                   | √ |
| 1d     | 700.10560 Post top for pivot post - 1 piece                                      | √ |
| 1e     | 700.10500 Standard post - 2 pieces                                               | √ |
| 1f     | 200.10720 Round blue cap - 3 pieces                                              | √ |
| 1g     | 000.03280 Allen key 6mm - 1 piece                                                | √ |
| 1h     | 700.10600 Clip-on level - 1 piece                                                | √ |
| 1i     | 700.10400 Top plate - 3 pieces                                                   | √ |
| 1j     | 700.10350 Round foot plate - 3 pieces                                            | √ |
| 1k     | 700.10230 Blue plastic strip - 6 pieces                                          | √ |
| 1l     | 200.11050 Flat blue cap - 1 piece                                                | √ |
| 1m     | 200.10780 Protector - 8 pieces                                                   | √ |
|        | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece              | o |

| Optional Protective Bag kit |                                              |   |
|-----------------------------|----------------------------------------------|---|
|                             | 002.10450 Easytrack transport bag - 2 pieces | o |

| Accessories |                                   |   |
|-------------|-----------------------------------|---|
| 2a          | 700.08310 Reacher - 61 cm (24 in) | o |
| 2b          | 700.08320 Reacher - 91 cm (36 in) | o |






### Floor to ceiling pressure fit system

## Easytrack 4 post with moving 3<sup>rd</sup> track

The Easytrack is a floor to ceiling pressure fit system requiring no permanent installation. Other versatile options are available such as 2 post, 3 post, bath tub adaptor and wall bracket configurations.



√ = Included / o = Optional

| # Product | Description        |
|-----------|--------------------|
| 9440030   | Easytrack - 4 post |

| # Part | Description                                                         |   |
|--------|---------------------------------------------------------------------|---|
| 1a     | 700.10250 X-Y Tracks 300 cm (118 in) - 3 pieces                     | √ |
| 1b     | 700.11005 Rail stopper - 2 pieces                                   | √ |
| 1c     | 700.10650 Standard trolley - 1 piece                                | √ |
| 1d     | 700.10701 Rail hook- 4 pieces                                       | √ |
| 1e     | 200.11050 Flat blue cap - 2 pieces                                  | √ |
| 1f     | 700.11270 X-Y mobile trolley - 1 piece                              | √ |
| 1g     | 200.10730 Drilling JIG - 1 piece                                    | √ |
| 1h     | 700.11280 X-Y fixed trolley - 1 piece                               | √ |
| 1i     | 700.10600 Clip-on level - 1 piece                                   | √ |
| 1j     | 000.03280 Allen key 6mm - 1 piece                                   | √ |
| 1k     | 000.02505 Allen key 4mm - 1 piece                                   | √ |
| 1l     | 200.10720 Round blue cap - 4 pieces                                 | √ |
| 1m     | 700.10400 Top plate - 4 pieces                                      | √ |
| 1n     | 700.10350 Round foot plate - 4 pieces                               | √ |
| 1o     | 700.10230 Blue plastic strip - 8 pieces                             | √ |
| 1p     | 700.10500 Standard post - 4 pieces                                  | √ |
|        | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece | o |

#### Optional Protective Bag kit

|                                                                    |                                    |   |
|--------------------------------------------------------------------|------------------------------------|---|
| 002.10450                                                          | Easytrack transport bag - 3 pieces | o |
| Custom bags are required if tracks are greater than 72 in / 183 cm |                                    |   |

#### Accessories

|    |                                   |   |
|----|-----------------------------------|---|
| 2a | 700.08310 Reacher - 61 cm (24 in) | o |
| 2b | 700.08320 Reacher - 91 cm (36 in) | o |



Floor to ceiling pressure fit system



Easytrack 1 post with a wall bracket


The Easytrack is a floor to ceiling pressure fit system requiring no permanent installation. Other versatile options are available such as 2 post, 3 post, 4 post, bath tub adaptor and wall bracket configurations.

| # Product | Description                            |
|-----------|----------------------------------------|
| 9211010   | Easytrack - 1 post with a wall bracket |

√ = Included / o = Optional

| # Part                      | Description                                                         |   |
|-----------------------------|---------------------------------------------------------------------|---|
| 1a                          | 700.10860 Wall Bracket - 100 cm (39 in) - 1 piece                   | √ |
| 1b                          | 700.10750 Extendable track & standard trolley - 1 piece             | √ |
| 1c                          | 700.10500 Standard post - 1 piece                                   | √ |
| 1d                          | 200.10720 Round blue cap - 1 piece                                  | √ |
| 1e                          | 700.10600 Clip-on level - 1 piece                                   | √ |
| 1f                          | 700.10400 Top plate - 1 piece                                       | √ |
| 1g                          | 700.10350 Round foot plate - 1 piece                                | √ |
| 1h                          | 700.10230 Blue plastic strip - 2 pieces                             | √ |
|                             | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece | o |
| Optional Protective Bag kit |                                                                     |   |
|                             | 002.10450 Easytrack transport bag - 1 piece                         | o |
| Accessories                 |                                                                     |   |
| 2a                          | 700.08310 Reacher - 61 cm (24 in)                                   | o |
| 2b                          | 700.08320 Reacher - 91 cm (36 in)                                   | o |







### Wall mounting system

## Easytrack 2 wall brackets system


Ceiling not suitable for a pressure fit mounting?  
Easytrack can also be fastened to the wall in various configurations.



| # Product | Description                 |
|-----------|-----------------------------|
| 9440030   | Easytrack - 2 wall brackets |

√ = Included / o = Optional

| # Part                      | Description                                                         |   |
|-----------------------------|---------------------------------------------------------------------|---|
| 1a                          | 700.10860 Wall Bracket - 100 cm (39 in) - 2 pieces                  | √ |
| 1b                          | 700.10750 Extendable track & standard trolley - 1 piece             | √ |
|                             | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece | o |
| Optional Protective Bag kit |                                                                     |   |
|                             | 002.10450 Easytrack transport bag - 1 piece                         | o |
| Accessories                 |                                                                     |   |
| 2a                          | 700.08310 Reacher - 61 cm (24 in)                                   | o |
| 2b                          | 700.08320 Reacher - 91 cm (36 in)                                   | o |




Floor to ceiling pressure fit system

Easytrack 2 post with a wall bracket

The Easytrack is a floor to ceiling pressure fit system requiring no permanent installation. Other versatile options are available such as 2 post, 3 post, 4 post and bath tub adaptor configurations.


| # Product | Description                            |
|-----------|----------------------------------------|
| 9321010   | Easytrack - 2 post with a wall bracket |



√ = Included / o = Optional

| # Part                      | Description                                                                      |   |
|-----------------------------|----------------------------------------------------------------------------------|---|
| 1a                          | 700.10860 Wall Bracket - 100 cm (39 in) - 1 piece                                | √ |
| 1b                          | 700.10850 Pivoting extendable track 300 cm (118 in) & standard trolley - 1 piece | √ |
| 1c                          | 700.10800 Extendable track 300 cm (118 in) & swivel trolley - 1 piece            | √ |
| 1d                          | 700.10500 Standard post - 2 pieces                                               | √ |
| 1e                          | 200.10720 Round blue cap - 2 pieces                                              | √ |
| 1f                          | 000.03280 Allen key 6mm - 1 piece                                                | √ |
| 1g                          | 700.10600 Clip-on level - 1 piece                                                | √ |
| 1h                          | 700.10400 Top plate - 2 pieces                                                   | √ |
| 1i                          | 700.10350 Round foot plate - 2 pieces                                            | √ |
| 1j                          | 700.10230 Blue plastic strip - 4 pieces                                          | √ |
| 1k                          | 200.11050 Flat blue cap - 1 piece                                                | √ |
|                             | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece              | o |
| Optional Protective Bag kit |                                                                                  |   |
|                             | 002.10450 Easytrack transport bag - 2 pieces                                     | o |
| Accessories                 |                                                                                  |   |
| 2a                          | 700.08310 Reacher - 61 cm (24 in)                                                | o |
| 2b                          | 700.08320 Reacher - 91 cm (36 in)                                                | o |






Wall mounting system

## Easytrack 3 wall brackets system with pivoting 2<sup>nd</sup> track


Ceiling not suitable for a pressure fit mounting?  
Easytrack can also be fastened to the wall in various configurations.



| # Product | Description                      |
|-----------|----------------------------------|
| 9303010   | Easytrack - Wall mounting system |

√ = Included / o = Optional

| # Part                      | Description                                                                      |   |
|-----------------------------|----------------------------------------------------------------------------------|---|
| 1a                          | 700.10860 Wall Bracket - 100 cm (39 in) - 3 pieces                               | √ |
| 1b                          | 700.10850 Pivoting extendable track 300 cm (118 in) & standard trolley - 1 piece | √ |
| 1c                          | 700.10800 Extendable track 300 cm (118 in) & swivel trolley - 1 piece            | √ |
| 1d                          | 000.03280 Allen key 6mm - 1 piece                                                | √ |
|                             | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece              | o |
| Optional Protective Bag kit |                                                                                  |   |
|                             | 002.10450 Easytrack transport bag - 2 pieces                                     | o |
| Accessories                 |                                                                                  |   |
| 2a                          | 700.08310 Reacher - 61 cm (24 in)                                                | o |
| 2b                          | 700.08320 Reacher - 91 cm (36 in)                                                | o |





Wall mounting system

Easytrack 4 wall brackets system  
with moving 3<sup>rd</sup> track

Ceiling not suitable for a pressure fit mounting?  
Easytrack can also be fastened to the wall in various configurations.


| # Product | Description                      |
|-----------|----------------------------------|
| 9404030   | Easytrack - Wall mounting system |

√ = Included / o = Optional



| # Part                      | Description                                                         |   |
|-----------------------------|---------------------------------------------------------------------|---|
| 1a                          | 700.10250 X-Y Tracks 300 cm (118 in) - 3 pieces                     | √ |
| 1b                          | 700.11005 Rail stopper - 2 pieces                                   | √ |
| 1c                          | 700.10650 Standard trolley - 1 piece                                | √ |
| 1d                          | 700.11270 X-Y mobile trolley - 1 piece                              | √ |
| 1e                          | 200.10730 Drilling JIG - 1 piece                                    | √ |
| 1f                          | 700.10600 Clip-on level - 1 piece                                   | √ |
| 1g                          | 000.03280 Allen key 6mm - 1 piece                                   | √ |
| 1h                          | 700.10701 Rail hook- 4 pieces                                       | √ |
| 1i                          | 700.10860 Wall Bracket - 100 cm (39 in) - 4 piece                   | √ |
|                             | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece | o |
| Optional Protective Bag kit |                                                                     |   |
|                             | 002.10450 Easytrack transport bag - 2 pieces                        | o |
| Accessories                 |                                                                     |   |
| 2a                          | 700.08310 Reacher - 61 cm (24 in)                                   | o |
| 2b                          | 700.08320 Reacher - 91 cm (36 in)                                   | o |







### Floor to ceiling pressure fit system

## Easytrack 2 post system with a bath adaptor

The Easytrack is a floor to ceiling pressure fit system requiring no permanent installation.

Other versatile options are available such as 2 post, 3 post, 4 post and wall bracket configurations.



| # Product | Description                            |
|-----------|----------------------------------------|
| 9220110   | Easytrack - 2 post with a bath adaptor |

√ = Included / o = Optional

| # Part                      | Description                                                         |   |
|-----------------------------|---------------------------------------------------------------------|---|
| 1a                          | 700.10880 Bath adaptor - 1 piece                                    | √ |
| 1b                          | 700.10750 Extendable track & standard trolley - 1 piece             | √ |
| 1c                          | 700.10455 Pivot post - 1 piece                                      | √ |
| 1d                          | 700.10500 Standard post - 1 piece                                   | √ |
| 1e                          | 200.10720 Round blue cap - 2 pieces                                 | √ |
| 1f                          | 700.10600 Clip-on level - 1 piece                                   | √ |
| 1g                          | 700.10400 Top plate - 2 pieces                                      | √ |
| 1h                          | 700.10350 Round foot plate - 1 pieces                               | √ |
| 1i                          | 700.10230 Blue plastic strip - 4 pieces                             | √ |
|                             | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece | o |
| Optional Protective Bag kit |                                                                     |   |
|                             | 002.10450 Easytrack transport bag - 2 pieces                        | o |
| Accessories                 |                                                                     |   |
| 2a                          | 700.08310 Reacher - 61 cm (24 in)                                   | o |
| 2b                          | 700.08320 Reacher - 91 cm (36 in)                                   | o |




Floor to ceiling pressure fit system

Easytrack 3 post system with a bath adaptor

The Easytrack is a floor to ceiling pressure fit system requiring no permanent installation. Other versatile options are available such as 2 post, 3 post, 4 post and wall bracket configurations.

| # Product | Description                            |
|-----------|----------------------------------------|
| 9330110   | Easytrack - 3 post with a bath adaptor |

√ = Included / o = Optional



| # Part                      | Description                                                         |   |
|-----------------------------|---------------------------------------------------------------------|---|
| 1a                          | 700.10880 Bath adaptor - 1 piece                                    | √ |
| 1b                          | 700.10750 Extendable track & standard trolley - 1 piece             | √ |
| 1c                          | 700.10250 X-Y Tracks 300 cm (118 in) - 1 piece                      | √ |
| 1d                          | 700.10500 Standard post - 2 pieces                                  | √ |
| 1e                          | 700.10455 Pivot post - 1 piece                                      | √ |
| 1f                          | 200.10720 Round blue cap - 3 pieces                                 | √ |
| 1g                          | 000.03280 Allen key 6mm - 1 piece                                   | √ |
| 1h                          | 700.10600 Clip-on level - 1 piece                                   | √ |
| 1i                          | 700.10400 Top plate - 3 pieces                                      | √ |
| 1j                          | 700.10350 Round foot plate - 2 pieces                               | √ |
| 1k                          | 700.10230 Blue plastic strip - 6 pieces                             | √ |
|                             | 700.10720 Extendable track (1.3-1.9m / 52-75 ") & trolley - 1 piece | o |
| Optional Protective Bag kit |                                                                     |   |
|                             | 002.10450 Easytrack transport bag - 3 pieces                        | o |
| Accessories                 |                                                                     |   |
| 2a                          | 700.08310 Reacher - 61 cm (24 in)                                   | o |
| 2b                          | 700.08320 Reacher - 91 cm (36 in)                                   | o |