

Legend Series L-325 Lift

Standard Features

- ADA Compliant Lift
- Stainless Steel Construction
- Powder coat finish
- Lift Capacity 325 lbs
- Battery Powered 24 volt system
- User Friendly Controls
- 5 Year Structural Warranty
- Arm rests
- Foot Rest
- Safety Seat Belt
- UPS Packaging for easy shipping
- Fits most standard pool applications
- UV resistant
- Right Seat Application
- Includes Battery Charging Station
- Key Lock out System
- Optional Drop in or Four Point Anchoring System



Global Lift Corp. Products are Laboratory tested to meet ADA Standards

Options



Protective Cover



Extra Battery



2014

New ADA Laws

- All of our lifts are compliant with the New ADA Laws
- All swimming pools in the USA that fall under Title II or Title III meet to comply with the new laws.

Title II (Public Industry) municipal pools, school pools, government owned pools, etc.


Title III (Private Industry) place of recreation, place of lodging, which addresses public accommodations

- Existing pools shall become compliant with accessibility by January 31, 2013
- There are tax incentives available
See details at www.GlobalLiftCorp.com

Legend
L-325

Deck Configurator Form

B



LEGEND MIN/MAX

A=(DRAFT) POOL DECK TO WATER LINE 0-7"

B=HEIGHT OF CURB 0-1"

Contact Global Lift Corp to confirm lift will work for your pool application.

Dealer: