

## Easy Lift Patient Lifting System - Bariatric

- Heavy-gauge steel construction
- Six-point spreader bar with 360° rotation tilts for extra comfort and safety
- Uses Lumex six-point, four-point slings or 2-point slings
- Unique, ergonomic foot pedal leg opening
- Large, easy to grip handles for easy maneuvering
- Lower base height maximizes stability and fits under most beds
- Rubber base coating protects furniture and walls
- 24V DC motor allows quiet smooth operation
- Ergonomically designed hand pendant with pendant hook
- Warning beep for low battery charge and weight capacity levels
- Easy access emergency stop button
- Emergency manual lowering device
- Available with maximum weight capacity of 600 lb
- White powder-coated finish
- Accessory storage bag included for slings and charger
- Optional digital scale available
- Designed to meet the requirements of HCPCS code:E0635
- 1 year limited warranty.

LF1090 1 ea



### Specifications

Type .....	Electric
Max. Weight Capacity.....	600 lb
Floor Clearance.....	5.75"
Lifting Range Min. Height.....	27.0"
Max. Height .....	71.0"
Overall Length.....	50.0"
Overall Width.....	26.2"
Max. Base Width .....	43.0"
Min. Base Width .....	26.2"
Weight Motor .....	137.0 lb
Warranty.....	24v DC
Frame .....	1 yr.
Battery.....	1 yr.

