

## Lumex® Easy Lift STS - LF2020

- Heavy-gauge steel construction.
- Designed to provide quick, safe assistance to those who have trouble standing.
- Unique, ergonomic foot pedal base opening.
- Large, easy to grip handles for easy maneuvering.
- Lower base height maximizes stability and fits under most beds.
- Rubber base coating protects furniture and walls.
- Adjustable knee pad unit.
- Quick release removable foot platform.
- Optional buttock strap for patients requiring additional support.
- 24V DC motor allows quiet, smooth operation.
- Ergonomically designed hand pendant with pendant hook.
- Warning beep for low battery charge and weight capacity levels.
- Easy access emergency stop button.
- Maximum weight capacity 400 lb.
- White powder-coated finish.
- Accessory storage bag included for slings and charger.
- Designed to meet the requirements of HCPCS code: E0635.
- 1 year limited warranty.

LF2020

1 ea



ITEM	Type	Max. Weight Capacity	Floor Clearance	Lifting Range Min. Height	Max. Height	Overall Length	Overall Width	Max. Base Width	Min. Base Width	Weight	Motor	Warranty Frame	Battery
LF2020	Electric Sit-Stand	400 lb	4.0"	33.0"	68.0"	40.0"	26.5"	37.0"	26.5"	116.6 lb	24v DC	1 yr.	1 yr.



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**www.grahamfield.com**