



redefining patient handling

Hoyer® Classics Product Overview

For more than 50 years, Hoyer lifters have defined the lifter category. From Ted Hoyer's early designs in the 1950's to more efficient models built with modern designs and materials, the Hoyer name has always stood for innovation. Hoyer's full line of truly functional and versatile lifts offers outstanding performance and style beyond compare!



HPL600WB/HPL600WBSC

- Power operated base
- Base clearance of 4.5"
- 4-point cradle design minimizes patient discomfort
- Long, padded handles offer various grip choices
- Floor pick-up capability
- Emergency stop, power manual lowering
- 600 lb. patient capacity



HPL402

- Adjustable U-base
- Reverse-mounted actuator for additional knee clearance
- 6-point non-sway cradle offers maximum sling options
- Built-in emergency stop button
- Manual release for emergency lowering
- Lift leg guards and caster bumpers
- 400 lb. patient capacity



P-C-HLA-2

- Powered lifting provides more focus on the patient being lifted
- Adjustable U-base with 5.5" clearance
- 6-point cradle works with an array of Hoyer slings
- Easy cleaning, durable chrome finish
- 400 lb. patient capacity

You must use Hoyer® slings with Hoyer lifts.



C-HLA-2

- Floor pick-up capability
- Adjustable legs
- Cradle design minimizes sway
- 6-point non-sway cradle offers maximum sling options
- Durable chrome finish maintains appearance
- Mast and base separate for easy travel or storage
- Standard 4" casters
- 400 lb. patient capacity

Hoyer® Classics Slings



**One Piece
Commode Sling with
Positioning Strap**

- Durable, polyester with closed-cell foam padding for comfort
- Aperture for toileting without removing sling
- Easy patient position adjustment
- Size range from S-L
- Weight Capacity: 600 lbs.



**One Piece
Sling with
Positioning Strap**

- Durable, polyester with closed-cell foam padding for comfort
- Provides easy position adjustment of patient
- Size range from S-L
- Weight Capacity: 600 lbs.



**U-Sling, Padded with
Head Support**

- Durable, polyester with closed-cell foam padding for comfort
- Nylon head support
- Full range of sizes from S-XL
- Weight Capacity: 600 lbs.



**U-Sling,
Padded**

- Durable, polyester with closed-cell foam padding for comfort
- Full range of sizes from S-XL
- Weight Capacity: 600 lbs.

You must use Hoyer® slings with Hoyer lifts.

Hoyer Classic Lifts	
HPL600WB/HPL600WBSC 600 lb. patient capacity	Power operated base, Base clearance of 4.5", 4-point cradle design minimizes patient sway
HPL402 400 lb. patient capacity	Adjustable U-base, Reverse-mounted actuator for additional knee clearance, 6-point cradle works with an array of Hoyer slings
P-C-HLA-2 400 lb. patient capacity	Adjustable U-base with 5.5" clearance, 6-point cradle works with an array of Hoyer slings
C-HLA-2 400 lb. patient capacity	Adjustable legs, Cradle design minimizes sway, 6-point non-sway cradle offers extra stability

Hoyer Classic Slings	X-Large	Large	Medium	Small
One Piece Commode Sling with Positioning Strap	—	70056	70057	70058
One Piece Sling with Positioning Strap	—	70051	70052	70053
U-Sling, Padded with Head Support	70010	70011	70012	70013
U-Sling, Padded	70000	70001	70002	70003

**For product and pricing questions,
call our customer care team today!**

Post Acute
(Nursing Homes,
LTACs, Assisted Living)
800.826.0270

Acute
(Hospitals and
Distribution Partners)
888.221.8461

Government
(VA and GSA)
800.966.6662



Joerns Healthcare Inc.
5001 Joerns Drive, Stevens Point WI 54481
(P) 800.826.0270 | (F) 800.457.8827
www.joerns.com | email: info@joerns.com