PAL Series



Descriptions for other variations:

- PAL Hi/Lo A variation of the PAL useful for facilities that have both inground pools and above ground spas.
- PAL Spa Designed for above ground pools and spas up to 48" above the pool deck.

Since its introduction in 1997, the PAL Portable Aquatic Lift has become the industry standard for providing access to swimming pools. The completely portable PAL is the perfect accessory for any commercial aquatic facility.

- Meets or exceeds Americans with Disabilities Act Access Guidelines (ADA/ABA 2004)
- 300 pound lifting capacity
- Completely portable no deck fixtures and easy to move by one person
- · Conveniently stores away when not in use
- Powered by a 24-volt rechargeable battery
- User operable via waterproof hand control
- 240° rotation ensures a safe transfer area
- Made from powder-coated stainless steel and aluminum

Model No.	Description	Shipping Wt.	L	W	Н
200-0000	PAL	1,100 lbs.	40"	48"	60"
250-0000	PAL Hi/Lo	1,100 lbs.	40"	48"	60"
275-0000	PAL Spa	1,100 lbs.	40"	48"	60"

Includes battery, charger, battery console cover, waterproof control, footrest and seatbelt assembly

Accessories:

- 170-1000 Arm Rest Assembly
- 920-2000 Total Cover
- 500-1000 Spineboard Attachment
- 900-2000 Stability Vest
- 100-1400 Wireless Controls











