



## Traxx Mobility Systems

The Titan 500 can reduce caregiver injury, improve the quality of life and provide a high standard of care with comfort and dignity.



## Traxx Mobility Systems TITAN 500

### TITAN 500 ADVANTAGES:

- Easy to use.
- Battery powered and rechargeable.
- Can lift patient from floor level if needed.
- Easy to move and re-install in new location.
- Installs in minutes.
- Lifts up to 500 lbs.
- No storage issues when not in use.

*Lifting the Standard of Home Care*

Traxx Mobility Systems

Customer Service

(855) TRAXXMS Or (855) 872-9967

9:00am to 5:00pm Eastern Time

Or visit the website at

[www.traxxms.com](http://www.traxxms.com)



*Lifting the Standard of Home Care*

33750 Freedom Road, Farmington, MI 48335 (248) 426-7700 Office (248) 426-7710 Fax

[www.traxxms.com](http://www.traxxms.com)

# Safe, reliable and affordable.

## IMPROVE THE QUALITY OF LIFE FOR BOTH THE LOVED ONE AND THE CAREGIVER.

The Titan 500 is an affordable alternative for compassionate home health care. Now you can have complete confidence when moving your loved one from bed to chair and back again.

### Lifting the Standard of Home Care

When it comes to individual care within the household, a bedroom is most often chosen as the primary care setting. Average room size is limited, usually less than 12' by 12' in total. With a bed and multiple pieces of furniture, the amount of space needed to properly maneuver the standard manual floor lift is extremely limited at best.

#### THE TRAXX TITAN 500

The Traxx Patient Lifting System improves the quality of life for individuals with limited mobility and their caregivers. The Traxx is safe, reliable, and extremely affordable. Due to its unique design, caregivers can lift and transfer their loved ones from bed to wheelchair and back again with ease, comfort, and security without the physical strains of manual lifting.

#### LIFT FEATURES

The Traxx Patient Lift is freestanding; there is no need to attach to walls or ceilings and installs in minutes. The lift motor is battery powered and easily recharged. The high strength aluminum frame can support up to 500 lbs and with the satin anodized finish makes it easy to clean. Lifting your loved one is as easy as one touch with the Traxx one-finger remote control.



#### Safe

The high strength aluminum frame can support up to 500 lbs, making it safe and secure.

#### Reliable

The ergonomic sling design combines strategically placed padding and adjustable strap attachment to a four point lift bar insuring a stable, secure and comfortable experience while being lifted or transferred.

#### Affordable

By dealing directly with the end user, Traxx is able to eliminate normal mark-ups and deliver to you an affordable and safe lift.



**TRAXX  
MOBILITY  
SYSTEMS**

Made In America  
Veteran Owned



CALL TOLL-FREE (855) TRAXXMS or (855) 872-9967

9:00am to 5:00pm Eastern Time and visit us on-line 24 hours a day

[www.traxxms.com](http://www.traxxms.com)

33750 Freedom Road, Ste. A, Farmington, MI 48335